

Certificate of Analysis

Product Information:

Orion Number 841108
Description 1000 PPM SODIUM STD 475 ML
Lot Number IZ1
Expiration Date 2 Years after Receipt of Product
Method of Analysis Potentiometric
Acceptance Range 1000 +/- 5 ppm Sodium as Sodium Chloride

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
CommandButton	1	1001 ppm Sodium as Sodium Chloride

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

Certificate Date: 01/19/04

Approved By:

Edward Henning

Quality Control Supervisor