

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	977011
<b>Description</b>	! ACID REAGENT 475 ML
<b>Lot Number</b>	IY1
<b>Expiration Date</b>	2 Years after Receipt of Product
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	4.2 +/- pH units

## Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	4.3	pH

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186 II	g
Potassium Dihydrogen Phosphate	186 I	g
Potassium Hydrogen Phthalate	185	h

**Certificate Date:** 02/17/04

**Approved By:**

*Edward Henning*

**Quality Control Supervisor**