

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	654204
<b>Description</b>	NON-IONIC SURF.TITRANT,475 ML
<b>Lot Number</b>	IZ1
<b>Expiration Date</b>	2 Years after Receipt of Product
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	0.050 M +/- 0.002 M BaCL2

## Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.051	Moles/Liter Barium Chloride

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

**Certificate Date:** 01/30/04

**Approved By:**

*Edward Henning*

**Quality Control Supervisor**