

Certificate of Analysis

Product Information:

Orion Number CAC45FZ000
Description 4000 NTU FORMAZIN STANDARD, 500ML
Lot Number IW1
Expiration Date 04-26-05
Method of Analysis Spectrophotometer

Actual Analysis:

<u>Analysis</u>	<u>Specification</u>	<u>Results</u>	<u>Unit of Measure</u>
1	4000+/-80 NTU FORMAZIN POLYMER	4080	ISO - NEPH METHOD

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
NOT APPLICABLE	NA	NA

Certificate Date: 05/03/04

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

245968-001 Rev. B

Thermo
ELECTRON CORPORATION