

Certificate of Analysis

Product Information:

Orion Number	013016A
Description	COND 2 CELL,SSTEEL K=0.1, 8 PIN F W/FLOW CELL
Lot Number	IY1018
Expiration Date	N/A
Method of Analysis	Conductivity
Acceptance Range	0.100 +/- 0.010 Cell Constant

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.094	Cell Constant (100uS Solution 011008 Lot: HP1)

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Chloride	2202	NA

Certificate Date: 05/10/04

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. B

Thermo
ELECTRON CORPORATION