

Certificate of Analysis

Product Information:

| | |
|---------------------------|-----------------------------------|
| Orion Number | 940909 |
| Description | TISAB II (ONE GALLON) |
| Lot Number | IT2 |
| Expiration Date | 2 Years after Receipt of Product |
| Method of Analysis | Potentiometric |
| Acceptance Range | pH 5.0 to 5.5/Fluoride < 0.10 ppm |

Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|------------------------|
| 1 | 5.46 | pH |
| 2 | 0.018 | ppm Fluoride |

Traceability:

| <u>N.I.S.T. Chemical</u> | <u>SRM #</u> | <u>Sublot</u> |
|--------------------------------|--------------|---------------|
| Disodium Hydrogen Phosphate | 186 II | f |
| Potassium Dihydrogen Phosphate | 186 I | f |
| Potassium Hydrogen Phthalate | 185 | g |
| Potassium Fluoride | 2203 | NA |

Certificate Date: 07/28/04

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

Thermo
ELECTRON CORPORATION