

Certificate of Analysis

Product Information:

Orion Number	013005A
Description	COND 4 CELL K=0.475, 8 PIN F
Lot Number	JU9006
Expiration Date	06/07
Method of Analysis	Conductivity
Acceptance Range	0.451 TO 0.498 Cell Constant

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.473	Cell Constant

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
NOT APPLICABLE	NA	NA

Certificate Date: 06/23/05

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. C

Thermo
ELECTRON CORPORATION