

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	013005A
<b>Description</b>	COND 4 CELL K=0.475, 8 PIN F
<b>Lot Number</b>	JU9011
<b>Expiration Date</b>	06/07
<b>Method of Analysis</b>	Conductivity
<b>Acceptance Range</b>	0.451 TO 0.498 Cell Constant

## Actual Analysis:


<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.467	Cell Constant

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
NOT APPLICABLE	NA	NA

**Certificate Date:** 06/23/05

**Approved By:**



**Quality Control Supervisor**

**Environmental Instruments**  
Water Analysis Instruments

166 Cummings Center

Beverly, MA  
01915

(978) 232-6000  
(978) 232-6015 fax

[www.thermo.com/water](http://www.thermo.com/water)  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. C

**Thermo**  
ELECTRON CORPORATION