

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	700902
<b>Description</b>	PURE WATER PH BUFFER C PH 9.15. 4 x 475 ML
<b>Lot Number</b>	JT1
<b>Expiration Date</b>	7/06
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	pH 9.19 +/- 0.05 pH Units @ 25.0 C

## Actual Analysis:


<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	9.20	pH @ 25 C

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186-II-g	NA
Potassium Dihydrogen Phosphate	186-I-g	NA
Sodium Tetraborate Decahydrate	187	d

**Certificate Date:** 07/21/05

**Approved By:**



**Quality Control Supervisor**

**Environmental Instruments**  
Water Analysis Instruments

166 Cummings Center

Beverly, MA  
01915

(978) 232-6000  
(978) 232-6015 fax

[www.thermo.com/water](http://www.thermo.com/water)  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. C

**Thermo**  
ELECTRON CORPORATION