

Certificate of Analysis

Product Information:

Orion Number	910001
Description	PH ELEC STORAGE SOLN, 475 ML
Lot Number	KZ1
Expiration Date	1/08
Method of Analysis	Potentiometric
Acceptance Range	3.00 M +/- 0.10 M Potassium Chloride, 6.75 +/- 0.2 pH units]

Actual Analysis:


<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	3.01	Moles/Liter Potassium Chloride
2	6.71	pH

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA
Disodium Hydrogen Phosphate	186-II-g	NA
Potassium Dihydrogen Phosphate	186 I	g
Sodium Bicarbonate	191	C
Sodium Carbonate	191	C

Certificate Date: 01/26/06

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. C

Thermo
ELECTRON CORPORATION