

Certificate of Analysis

Product Information:

Orion Number 700402
Description PURE WATER PH BUFFER B PH 4.10, 4 X 475 ML
Lot Number KR1
Expiration Date 09/08
Method of Analysis Potentiometric
Acceptance Range 0.010 M +/- 0.001 M Potassium Chloride / pH 4.10 +/- 0.03 pH units

Actual Analysis:


<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.010	Moles/Liter Potassium Chloride
2	4.13	pH

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA
Disodium Hydrogen Phosphate	186 II	g
Potassium Dihydrogen Phosphate	186 I	g
Potassium Hydrogen Phthalate	185	h

Certificate Date: 09/15/06

Approved By:



Quality Representative

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. D

Thermo
ELECTRON CORPORATION