

Certificate of Analysis

Product Information:

Orion Number	941609
Description	! SAOB REAGENT PACK OF 4 BOTTLES FOR SULFIDE
Lot Number	KR1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	Std. A 2.0 +/- 0.5 M NaOH / 0.20 +/- 0.05 M Na2EDTA Std. B 0.100 +/- 0.015 M Ascorbic Acid

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	1.74	Moles/Liter Sodium Hydroxide
2	0.095	Moles/Liter Ascorbic Acid

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Hydrogen Phthalate	185H	NA
Sodium Chloride	2201	NA

Certificate Date: 09/27/06

Approved By:



Quality Representative

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

Thermo
ELECTRON CORPORATION