

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	700702
<b>Description</b>	PURE WATER PH BUFFER A PH 6.97, 4 X 475 ML
<b>Lot Number</b>	LP1
<b>Expiration Date</b>	11/09
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	0.0100 +/- 0.0005 M Potassium Chloride/ 6.97 +/- 0.03 pH Units

## Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	.0101	Moles/Liter Potassium Chloride

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

**Certificate Date:** 11/28/07

**Approved By:**



**Quality Representative**

**Environmental Instruments**  
Water Analysis Instruments

166 Cummings Center

Beverly, MA  
01915

(978) 232-6000  
(978) 232-6015 fax

www.thermo.com/water  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. D

**Thermo**  
ELECTRON CORPORATION