



Conductivity Values at Temperature for Orion  
111.9 mS/cm Conductivity Standard

Temperature (degrees C)	Conductivity Standard Orion 011005 (mS/cm)
0	65.10
1	66.84
2	68.59
3	70.35
4	72.12
5	73.91
6	75.70
7	77.50
8	79.32
9	81.15
10	82.98
11	84.83
12	86.69
13	88.56
14	90.45
15	92.34
16	94.24
17	96.15
18	98.08
19	100.0
20	102.0
21	103.9
22	105.9
23	107.9
24	109.9
25	111.9
26	113.9
27	115.9
28	117.9
29	120.0
30	122.0
31	124.1
32	126.2
33	128.3
34	130.4
35	132.5
36	134.6
37	136.7
38	138.9
39	141.0
40	143.2
41	145.4
42	147.6
43	149.8
44	152.0
45	154.2
46	156.4
47	158.7
48	160.9
49	163.2
50	165.4