

Certificate of Analysis

Product Information:

Orion Number	510011
Description	NO CAL PH ELEC FILL SOLN 5 X 60ML
Lot Number	MW1
Expiration Date	04/2010
Method of Analysis	Potentiometric
Acceptance Range	PH.7.050 +/-0.015 PH units

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	7.055	pH @ 25 degrees C
2	4.00	Moles/Liter Potassium Iodide
3	0.0070	Moles/Liter Iodine

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186-II-g	NA
Potassium Dihydrogen Phosphate	186 I	g
Sodium Bicarbonate	191-I-c	NA
Sodium Bicarbonate	191	C
Sodium Chloride	2201	NA
Potassium Dichromate	136e	NA

Certificate Date: 05/02/08

Approved By:



Quality Representative

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. D

Thermo
ELECTRON CORPORATION