

Certificate of Analysis

Product Information:

Orion Number 951211
Description ! AMMONIA PH ADJST ISA 475ML
Lot Number MW1
Expiration Date N/A; Recommendation usage 2 years from date of receipt
Method of Analysis Potentiometric
Acceptance Range 5.00+/-0.25 M Sodium Hydroxide/ 0.05 M +/-0.05 Disodium EDTA

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	4.95	Moles/Liter Sodium Hydroxide
2	0.05	Moles/Liter Disodium E.D.T.A.

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA
Calcium Carbonate	915b	NA

Certificate Date: 05/06/08

Approved By:



Quality Representative

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

Thermo
ELECTRON CORPORATION