

Certificate of Analysis

Product Information:

| | |
|---------------------------|---|
| Orion Number | 954606 |
| Description | Standard Soln |
| Lot Number | MU1 |
| Expiration Date | 06/2009 |
| Method of Analysis | Potentiometric |
| Acceptance Range | 0.1000 M +/- 0.0008 M NaNO ₂ |

Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|----------------------------|
| 1 | 0.1003 | Moles/Liter Sodium Nitrite |

Traceability:

| <u>N.I.S.T. Chemical</u> | <u>SRM #</u> | <u>Sublot</u> |
|--|--------------|---------------|
| TitriStar Ceric Sulfate 0.1 N "tested against NIST SRM Arsenic Trioxide" | | 83d |

Certificate Date: 07/02/08

Approved By:



Quality Representative

Environmental Instruments
Water Analysis Instruments

166 Cummings Center

Beverly, MA
01915

(978) 232-6000
(978) 232-6015 fax

www.thermo.com/water
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

244721-001 Rev. D

Thermo
ELECTRON CORPORATION