

# Shandon Embedding Accessories

## Shandon PARA TRIMMER®



The Shandon PARA TRIMMER is designed to remove excess wax from the outside of a tissue cassette after embedding by sliding the outside of the block across the heated surface.

Paraffin should be removed to properly seat the block in the quick-release clamp or vise to achieve the best section quality from your microtome.

- Heated surface with ridges.
- Multiple blocks may be passed across the surface at one time.
- Patented, compact, dual voltage unit.

Description	Order Number
<b>Shandon PARA TRIMMER</b> 110-240V, 50/60Hz	B3120205

### Physical Specifications

#### Dimensions (H x W x D)

6.75" x 5.25" x 7.125" (17.2 x 18.5 x 13.3 cm)

#### Weight

2 lbs. (1 kg)

#### Power Requirements

110-220V 50/60Hz (120VA)

### Agency Approvals

Designed to meet IEC 61010, CSA, UL, and CE Mark.

## Shandon Paraffin



Shandon paraffins are available in a wide range of formulations to meet virtually any customer need.

- Paraplast®
- Paraplast Plus®
- Paraplast X-tra®
- Shandon Precision Cut
- Shandon Histoplast

\* Refer to Anatomical Pathology Product Catalog page 11 for complete details.

## Shandon Wax Dispenser



- Holds up to 5 liters of molten wax at a thermostatically controlled temperature.
- Separate tap heater ensures an even flow of wax.
- Convenient as a stand-alone or auxiliary source of molten wax.
- Stainless steel outer casing and lid.
- Aluminum inner reservoir for good thermal conductivity and improved temperature control.
- Non-drip lever tap controls wax delivery and can be locked in the open position if required.
- Reservoir filter screen prevents coarse particles from blocking the delivery tap.

Description	Order Number
<b>Shandon Wax Dispenser</b> 110-120V, 50/60Hz	3120069
220-240V, 50/60Hz	3120057

### Physical Specifications

#### Dimensions (H x W x D)

Diameter: 8.46" x Height: 16.1" (21.5 x 41 cm)

#### Weight

8 lbs. (3.7 kg)

#### Power Requirements

110-120V or 220-240V, 50/60 Hz, 250 Watts AC

#### Temperature Range

40 - 75°C, ±2°C

### Included with Product\*

Operator Manual 312069

\* Contents, Packaging, and Quantities of items listed under "Included with Product" are subject to change without prior notice.

### Agency Approvals

Designed to meet IEC 61010, CSA, UL, and CE Mark.

## Shandon Paraffin Section Flotation Bath

## Shandon Circular Paraffin Section Flotation Bath



- Holds water at a constant temperature to flatten wax or resin sections prior to lifting them onto a microscope slide.
- Also used for hot staining except with corrosive substances. The top rim of the unit has generous space for drying slides.
- Controls have been recessed to provide the user with an added measure of safety.
- A drop-on stainless steel lid is provided to cover the bowl when not in use.
- Deep tapered bowl is made from stainless steel to offer maximum resistance against corrosion.
- Both the bowl and the upper casing have a matte-black PTFE coating providing visual contrast against wax sections.
- Lower steel casing has a resistant epoxy coating.
- Power inlet socket at the back of the unit has fuses on both the negative and positive lines for optimum protection.

- Compact unit keeps water at a constant temperature to flatten wax or resin sections prior to lifting them onto a microscope slide.
- It may also be used for hot staining except with corrosive substances.
- The bath is made from aluminum with an anodized matte-black finish. It has a wide rim with ample space for slide drying.
- The power cable has molded connections at both terminals.
- The power inlet socket at the back of the unit has fuses on both the negative and positive lines for optimum protection.

Description	Order Number
<b>Shandon Paraffin Section Flotation Bath</b>	110-120V 50/60Hz 3120070
	220-240V 50/60Hz 3120058

### Physical Specifications

#### Dimensions (H x W x D)

External: 4.2" x 10.8" x 12.8" (10.7 x 27.5 x 32.5 cm)

Bowl: 8.7" Dia. x 3.34" Deep (22.2 x 8.5 cm)

#### Weight

8.26 lbs. (3.75 kg)

#### Temperature Range

30° - 65°C (± 1.5°) at 20°C ambient

#### Power Requirements

110-120V or 220-240V, 50/60Hz, 250 VA

Description	Order Number
<b>Shandon Circular Paraffin Section Flotation Bath</b>	110-120V 50/60Hz 3120071
	220-240V 50/60Hz 3120059

### Physical Specifications

#### Dimensions (H x W x D)

External: 12.5" Dia. x 6.26" Deep (32 x 15.9 cm)

Bowl: 8.66" Dia. x 2.95" Deep (22 x 7.5 cm)

#### Weight

5.5 lbs. (2.5 kg)

#### Temperature Range

30° - 70°C (± 1.5°) at 20°C ambient

#### Power Requirements

110-120V or 220-240V, 50/60Hz, 250 VA

Included with Product*	Order Number
Operator Manual	312071

\* Contents, Packaging, and Quantities of items listed under "Included with Product" are subject to change without prior notice.

### Agency Approvals

CE Mark.

Included with Product*	Order Number
Operator Manual	312070

\* Contents, Packaging, and Quantities of items listed under "Included with Product" are subject to change without prior notice.

### Agency Approvals

CE Mark.

**Anatomical Pathology USA**  
Tel: 1-800-547-7429 1-412-788-1133  
Fax: 1-412-788-1138

**Anatomical Pathology International**  
Tel: 0800 0189396 +44 (0) 1928 562600  
Fax: +44 (0) 1928 562627

**Anatomical Pathology Latin America/Canada**  
Tel: 1-412-788-2460  
Fax: 1-412-788-2480

Email: [clinical.diagnostics@thermo.com](mailto:clinical.diagnostics@thermo.com)

Internet: [www.thermo.com/shandon](http://www.thermo.com/shandon)