

NEW! Sorvall® RC-6™ Plus Superspeed Centrifuge

WATCH your  
productivity soar!

The Sorvall RC-6 Plus from Thermo is a powerful, compact centrifuge that combines superior functionality, unparalleled safety and proven reliability in a user-friendly design. It's built to ensure consistent, high-quality separations for your high-volume processing requirements.

- **High Speed, High Capacity Processing**

Offering 4 liters of capacity and processing speeds of up to 22,000 rpm (55,200 x g), the RC-6 Plus enables extremely fast processing of your samples.

- **Maximum Sample Protection**

The RC-6 Plus spins samples at atmospheric pressure, so there is less risk of tube or bottle leakage and no vacuum seals or filters to maintain.

- **Safer by Design**

The RC-6 Plus meets CE and CSA safety standards without the need to anchor the unit to the floor.

- **Versatile**

Choose from a range of carbon fiber, aluminum and continuous-flow rotors to meet your RCF, application and sample volume requirements.



**SPECIAL  
PACKAGE PRICING!**

**Save up to \$15,000**  
on a Sorvall RC-6 Plus-  
FIBERLite® rotor package!

*More details on reverse side.*

# DISCOVER the SORVALL *difference*

## Safety, reliability and performance – an unbeatable combination

The new Sorvall RC-6 Plus superspeed centrifuge is the **perfect choice for high-volume sample processing**, offering HIGH CAPACITY, HIGH SPEED, and HIGH G-FORCES along with advanced functionality that will **keep your research on the cutting edge.**



### NEW! FIBERLite Carbon Fiber Rotors –

High-performance carbon fiber rotors from FIBERLite provide a lighter, safer, longer-lasting alternative to conventional metallic rotors. FIBERLite adapters let you spin centrifuge tubes up to maximum speeds – reducing processing time and saving steps.

### Unrivaled Functionality –

Intuitive controls on the RC-6 Plus provide access to powerful functionality, such as step-run,  $v^2/dt$  integrator, automatic RCF/RPM conversion, start/stop time, and custom programs. The user lockout feature limits centrifuge access to authorized users in a multi-user environment.

### Unparalleled Safety –

The RC-6 Plus meets CE and CSA safety requirements without the need to anchor the instrument to the floor. And, by spinning samples at atmospheric pressure—without the need for a vacuum—the RC-6 Plus minimizes maintenance and helps prevent sample leakage, rotor imbalance and run shut-downs.

## Save up to \$15,000!

Purchase a Sorvall RC-6 Plus-FIBERLite rotor package – and SAVE. Package includes the RC-6 Plus and your choice of TWO FIBERLite carbon-fiber rotors!

FIBERLite Fixed Angle Rotors	Capacity (place x mL)	Max. Speed (rpm)	Max RCF (x g)	Cat. No.
SLC-4000	4 x 1,000	8,000	12,006	76980
SLC-3000	6 x 500	11,000	20,449	78510
F14S	6 x 250	14,000	25,862	78500
F21S	8 x 50	18,000	38,397	46923
F13S	14 x 50	12,500	26,727	46922

Contact your Thermo representative at  
**1-866-9-THERMO (984-3766)**  
or visit us online at  
**[www.thermo.com/centrifuge](http://www.thermo.com/centrifuge)**

Offer valid in U.S. only. Offer expires 12/15/2006.  
Not valid with any other offers or agreements.  
Void where prohibited. Bonus offers based on U.S. list prices.