

Material Safety Data Sheet
Factor II Deficient Plasma
Section 1 - Chemical Product and Company Identification
MSDS Name:

Factor II Deficient Plasma

Catalog Numbers:

102000, 176404, 22176404, 252000

Synonyms:

None.

Company Identification:
 Fisher Diagnostics
 Fisher Scientific Company LLC
 8365 Valley Pike
 Middletown, VA 22645-0307
Company Phone Number:

(800) 528-0494

Emergency Phone Number:

(800) 528-0494

CHEMTREC Phone Number:

800-424-9300

CHEMTREC (International):

703-527-3887

Section 2 - Composition, Information on Ingredients

CAS#	Chemical Name:	Percent	EINECS/ ELINCS
Not available	Bovine Plasma Citrate	Balance	Not available
7727-43-7	Barium sulfate	<0.1%	231-784-4
Not available	Factor II Deficient Plasma	35.0%	Not available
68-04-2	Trisodium citrate anhydrous	3.8	200-675-3
7365-45-9	HEPES	1.2%	230-907-9
7447-41-8	Lithium chloride	0.4%	231-212-3
1332-58-7	Kaolin	0.3%	Not available
9087-70-1	Aprotinin, Bovine	0.3%	232-994-9
9011-04-5	Polybrene	0.2%	Not available
1310-73-2	Sodium hydroxide	0.1%	215-185-5

C

35

Material Safety Data Sheet
Factor II Deficient Plasma**Section 3 - Hazards Identification****EMERGENCY OVERVIEW**

Appearance: white solid

Caution! Each unit of source material used in the preparation of this product has been tested by an FDA licensed method and found non-reactive for HBsAg and negative for antibodies to HIV and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS, or other infectious diseases. This product, like all materials of human origin, should be handled as potentially infectious biological material. May cause irritation. May cause central nervous system effects.

Target Organs: central nervous system.

Potential Health Effects**Eye:**

May cause eye irritation.

Skin:

Prolonged and/or repeated contact may cause irritation and/or dermatitis. Product tests non-reactive for HBsAg and negative for antibodies to HIV and HCV. However, material should be handled as potentially infectious.

Ingestion:

May cause irritation of the digestive tract.

Inhalation:

Prolonged exposure may result in dizziness and general weakness.

Chronic:

No information found.

Section 4 - First Aid Measures**Eyes:**

Flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical aid immediately.

Skin:

Get medical aid. Flush skin with plenty of soap and water.

Ingestion:

If victim is conscious and alert, give 2-4 cupfuls of milk or water. Get medical aid.

Inhalation:

Get medical aid immediately. Remove from exposure and move to fresh air immediately. If not breathing, give artificial respiration. If breathing is difficult, give oxygen.

Notes to Physician:

Treat symptomatically and supportively.

Material Safety Data Sheet

Factor II Deficient Plasma

Section 5 - Fire Fighting Measures

General Information:

As in any fire, wear a self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or equivalent), and full protective gear.

Extinguishing Media:

For small fires, use water spray, dry chemical, carbon dioxide or chemical foam.

Autoignition Temperature:

Not applicable.

Explosion Limits:

Lower: Not available. Upper: Not available.

Flash Point:

Not applicable.

NFPA Rating:

(estimated) Health: 1; Flammability: 0; Reactivity: 0

Section 6 - Accidental Release Measures

General Information:

Use proper personal protective equipment as indicated in Section 8.

Spills/Leaks:

In case of spill, sweep up solid or absorb liquid with inert material, then place into a biohazard container. Wear protective equipment.

Section 7 - Handling and Storage

Handling:

Wash thoroughly after handling. Use with adequate ventilation. Avoid contact with skin and eyes. Avoid ingestion and inhalation.

Storage:

Store in a cool place in the original container and protect from sunlight.

Section 8 - Exposure Controls, Personal Protection

Engineering Controls:

Good general ventilation should be sufficient to control airborne levels.

Exposure Limits

Chemical Name:	ACGIH	NIOSH	OSHA
Barium sulfate	10 mg/m ³ TWA	10 mg/m ³ TWA (total dust); 5 mg/m ³ TWA (respirable dust)	15 mg/m ³ TWA (total dust); 5 mg/m ³ TWA (respirable fraction)
Trisodium citrate anhydrous	None of the components are on this list.	None of the components are on this list.	None of the components are on this list.

Material Safety Data Sheet

Factor II Deficient Plasma

HEPES	None of the components are on this list.	None of the components are on this list.	None of the components are on this list.
Lithium chloride	None of the components are on this list.	None of the components are on this list.	None of the components are on this list.
Kaolin	2 mg/m ³ TWA (respirable fraction, particulate matter containing no asbestos and	10 mg/m ³ TWA (total dust); 5 mg/m ³ TWA (respirable dust)	15 mg/m ³ TWA (total dust); 5 mg/m ³ TWA (respirable fraction)
Aprotinin, Bovine	None of the components are on this list.	None of the components are on this list.	None of the components are on this list.
Polybrene	None of the components are on this list.	None of the components are on this list.	None of the components are on this list.
Sodium hydroxide	2 mg/m ³ Ceiling	10 mg/m ³ IDLH 2 mg/m ³ Ceiling	2 mg/m ³ TWA

OSHA Vacated PELs

Barium sulfate: 10 mg/m³ TWA (total dust); 5 mg/m³ TWA (respirable fraction)

Kaolin: 10 mg/m³ TWA (total dust); 5 mg/m³ TWA (respirable fraction)

Personal Protective Equipment

Eyes:

Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

Skin:

Wear appropriate protective gloves to prevent skin exposure.

Clothing:

Wear appropriate protective clothing to prevent skin exposure.

Respirators:

Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN 149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.

Section 9 - Physical and Chemical Properties

Physical State: Solid

Color: white

Odor: Odorless.

pH: No information found

Vapor Pressure: No information found

Vapor Density: No information found

Evaporation Rate: No information found

Viscosity: No information found

Boiling Point: Not available.

Freezing/Melting Point: No information found

Decomposition Temperature: No information found

Solubility in water: Not available.

Specific Gravity/Density: No information found

Molecular Formula: Not applicable.

Molecular Weight: No information found

Material Safety Data Sheet
Factor II Deficient Plasma**Section 10 - Stability and Reactivity****Chemical Stability:**

Stable under normal temperatures and pressures.

Conditions to Avoid:

None reported.

Incompatibilities with Other Materials

None reported.

Hazardous Decomposition Products

None.

Hazardous Polymerization

Has not been reported

Section 11 - Toxicological Information**RTECS:**

CAS# 7727-43-7: CR0600000; CR0630000.

CAS# 68-04-2: GE8300000.

CAS# 7365-45-9: TL6809000.

CAS# 7447-41-8: OJ5950000.

CAS# 1332-58-7: GF1670500.

CAS# 9087-70-1: YN5080000.

CAS# 9011-04-5: MO1400000.

CAS# 1310-73-2: WB4900000.

LD50/LC50:

CAS# 7727-43-7:

No information foundCAS# 68-04-2:

No information foundCAS# 7365-45-9:

No information foundCAS# 7447-41-8:

Draize test, rabbit, eye: 100 mg/24H Moderate

Draize test, rabbit, skin: 500 mg/24H Severe

Oral, mouse: LD50 = 1165 mg/kg

Oral, mouse: LD50 = 1165 mg/kg

Oral, rabbit: LD50 = 800 mg/kg

Oral, rabbit: LD50 = 775 mg/kg

Oral, rat: LD50 = 526 mg/kg

Oral, rat: LD50 = 1530 mg/kg

Skin, rabbit: LD50 = 1629 mg/kg

Skin, rat: LD50 = 1488 mg/kg.

CAS# 1332-58-7:

No information foundCAS# 9087-70-1:

No information foundCAS# 9011-04-5:

No information foundCAS# 1310-73-2:

Draize test, rabbit, eye: 400 ug Mild

Draize test, rabbit, eye: 1% Severe

Draize test, rabbit, eye: 50 ug/24H Severe

Draize test, rabbit, eye: 1 mg/24H Severe

Draize test, rabbit, skin: 500 mg/24H Severe.

Material Safety Data Sheet
Factor II Deficient Plasma**Carcinogenicity:**

CAS# 7727-43-7: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 68-04-2: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 7365-45-9: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 7447-41-8: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 1332-58-7: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 9087-70-1: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 9011-04-5: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.
CAS# 1310-73-2: Not listed as a carcinogen by ACGIH, IARC, NIOSH, NTP, OSHA, or CA Prop 65.

Epidemiology:

No information found

Teratogenicity:

No information found

Reproductive:

No information found

Mutagenicity:

No information found

Neurotoxicity:

No information found

Other:

See actual entry in RTECS for complete information.

Section 12 - Ecological Information

No information found

Section 13 - Disposal Considerations

Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. US EPA guidelines for the classification determination are listed in 40 CFR Part 261.3. Additionally, waste generators must consult state and local hazardous waste regulations to ensure complete and accurate classification.

RCRA P Series Wastes

None of the components are on this list.

RCRA U Series Wastes

None of the components are on this list.

Section 14 - Transport Information**US DOT****Canadian TDG**

Proper Shipping Name: No information found

No information found

Hazard Class:

UN Number:

Packing Group:

Material Safety Data Sheet
Factor II Deficient Plasma**Section 15 - Regulatory Information****US Federal****TSCA**

CAS# 7727-43-7 is listed on the TSCA Inventory.

CAS# 68-04-2 is listed on the TSCA Inventory.

CAS# 7365-45-9 is listed on the TSCA Inventory.

CAS# 7447-41-8 is listed on the TSCA Inventory.

CAS# 1332-58-7 is listed on the TSCA Inventory.

CAS# 9087-70-1 is not listed on the TSCA Inventory. It is for research and development use only.

CAS# 9011-04-5 is not listed on the TSCA Inventory. It is for research and development use only.

CAS# 1310-73-2 is listed on the TSCA Inventory.

Health and Safety Reporting List

None of the components are on this list.

Chemical Test Rules

None of the components are on this list.

TSCA Section 12b

None of the components are on this list.

TSCA Significant New Use Rule (SNUR)

None of the components are on this list.

CERCLA Hazardous Substances and corresponding RQs

CAS# 1310-73-2: 1000 lb final RQ; 454 kg final RQ

SARA Section 302 Extremely Hazardous Substances

None of the components are on this list.

SARA Hazard Categories

CAS# 7727-43-7: acute.

CAS# 7365-45-9: acute.

CAS# 7447-41-8: acute, chronic.

CAS# 1310-73-2: acute, reactive.

CERCLA/SARA Section 313

This chemical is not at a high enough concentration to be reportable under Section 313.

Clean Air Act - Hazardous Air Pollutants (HAPs)

None of the components are on this list.

Clean Air Act - Class 1 Ozone Depletors

None of the components are on this list.

Clean Air Act - Class 2 Ozone Depletors

None of the components are on this list.

Clean Water Act - Hazardous Substances

CAS# 1310-73-2 is listed as a Hazardous Substance under the CWA.

Clean Water Act - Priority Pollutants

None of the components are on this list.

Clean Water Act - Toxic Pollutants

None of the components are on this list.

Material Safety Data Sheet
Factor II Deficient Plasma**OSHA - Highly Hazardous**

None of the components are on this list.

US State**State Right to Know**

Barium sulfate can be found on the following state Right-to-Know lists: New Jersey, (listed as Barium compounds, n.o.s.), Pennsylvania, Minnesota, Massachusetts.

Kaolin can be found on the following state Right-to-Know lists: Pennsylvania, Minnesota, Massachusetts.

Sodium hydroxide can be found on the following state Right-to-Know lists: California, New Jersey, Pennsylvania, Minnesota, Massachusetts.

California Regulations

CAS# Not available: No information found
CAS# 7727-43-7: No information found
CAS# Not available: No information found
CAS# 68-04-2: No information found
CAS# 7365-45-9: No information found
CAS# 7447-41-8: No information found
CAS# 1332-58-7: No information found
CAS# 9087-70-1: No information found
CAS# 9011-04-5: No information found
CAS# 1310-73-2: No information found

European/International Regulations**European Labelling in Accordance with EC Directives:**

Hazard Symbols:

Risk Phrases:

Safety Phrases:

WGK (Water Danger/Protection)

CAS# Not available: No information found
CAS# 7727-43-7: 0
CAS# Not available: No information found
CAS# 68-04-2: 0
CAS# 7365-45-9: No information found
CAS# 7447-41-8: 1
CAS# 1332-58-7: 0
CAS# 9087-70-1: No information found
CAS# 9011-04-5: 2
CAS# 1310-73-2: 1

United Kingdom Occupational Exposure Limits

No information found

Canadian DSL/NDSL

CAS# 7727-43-7 is listed on Canada's DSL List.
CAS# 68-04-2 is listed on Canada's DSL List.
CAS# 7365-45-9 is listed on Canada's DSL List.
CAS# 7447-41-8 is listed on Canada's DSL List.
CAS# 1332-58-7 is listed on Canada's DSL List.
CAS# 1310-73-2 is listed on Canada's DSL List.

Canadian WHMIS Classifications

This product does not have a WHMIS classification.

Material Safety Data Sheet
Factor II Deficient Plasma**Canadian Ingredient Disclosure List**

CAS# Not available is not listed on Canada's Ingredient Disclosure List.
CAS# 7727-43-7 is not listed on Canada's Ingredient Disclosure List.
CAS# Not available is not listed on Canada's Ingredient Disclosure List.
CAS# 68-04-2 is not listed on Canada's Ingredient Disclosure List.
CAS# 7365-45-9 is not listed on Canada's Ingredient Disclosure List.
CAS# 7447-41-8 is not listed on Canada's Ingredient Disclosure List.
CAS# 1332-58-7 is not listed on Canada's Ingredient Disclosure List.
CAS# 9087-70-1 is not listed on Canada's Ingredient Disclosure List.
CAS# 9011-04-5 is not listed on Canada's Ingredient Disclosure List.
CAS# 1310-73-2 is listed on Canada's Ingredient Disclosure List.

Exposure Limits

CAS# 7727-43-7: OEL-AUSTRALIA:TWA 10 mg/m³
OEL-BELGIUM:TWA 10 mg/m³ (resp. dust)
OEL-UNITED KINGDOM:TWA 2 mg/m³
OEL IN BULGARIA, COLOMBIA, JORDAN, KOREA check ACGIH TLV
OEL IN NEW ZEALAND, SINGAPORE, VIETNAM check ACGI TLV
CAS# 1332-58-7: OEL-AUSTRALIA:TWA 10 mg/m³
OEL-BELGIUM:TWA 10 mg/m³
OEL-FRANCE:TWA 10 mg/m³
OEL-THE NETHERLANDS:TWA 10 mg/m³
OEL-SWITZERLAND:TWA 6 mg/m³
OEL IN BULGARIA, COLOMBIA, JORDAN, KOREA check ACGIH TLV
OEL IN NEW ZEALAND, SINGAPORE, VIETNAM check ACGI TLV
CAS# 1310-73-2: OEL-AUSTRALIA:TWA 2 mg/m³
OEL-BELGIUM:STEL 2 mg/m³
OEL-DENMARK:TWA 2 mg/m³
OEL-FINLAND:TWA 2 mg/m³
OEL-FRANCE:TWA 2 mg/m³
OEL-GERMANY:TWA 2 mg/m³
OEL-JAPAN:STEL 2 mg/m³
OEL-THE NETHERLANDS:TWA 2 mg/m³
OEL-THE PHILIPPINES:TWA 2 mg/m³
OEL-SWEDEN:TWA 2 mg/m³
OEL-SWITZERLAND:TWA 2 mg/m³;STEL 4 mg/m³
OEL-THAILAND:TWA 2 mg/m³
OEL-TURKEY:TWA 2 mg/m³
OEL-UNITED KINGDOM:TWA 2 mg/m³;STEL 2 mg/m³
OEL IN BULGARIA, COLOMBIA, JORDAN, KOREA check ACGIH TLV
OEL IN NEW ZEALAND, SINGAPORE, VIETNAM check ACGI TLV

Section 16 - Other Information

Color information has been
MSDS Creation Date: September 2, 1997
Revision Date: August 1, 2005



Material Safety Data Sheet
Factor II Deficient Plasma

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall the company be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential, or exemplary damages howsoever arising, even if the company has been advised of the possibility of such damages.