

Aspire RP30 Desalting Tips: Pipette Compatibility Chart

MANUFACTURER	MODEL	VOLUME (µL)
THERMO SCIENTIFIC FINNPIPETTE	FOCUS	10-100
	FOCUS	30-300
	NOVUS ELECTRONIC	30-300
BIOHIT	PROLINE	20-200
	PROLINE ELECTRONIC	5-100
	TRANSFERPETTE	20-200
BRANDTECH	SOFTLINE	10-100
	SOFTLINE	25-200
CAPPAERO	CAPPAERO	25-200
COSTAR	ONEPETTE	10-100
	ONEPETTE	20-200
EPPENDORF	PHYSIO CARE RESEARCH	10-100
	PHYSIO CARE RESEARCH	20-200
	PHYSIO CARE RESEARCH	50-200
	PHYSIO CARE RESEARCH	10-100
	2000 RESEARCH	50-200
	2100 RESEARCH	20-200
FISHER BRAND	FINNPIPETTE	40-200
	FINNPIPETTE II	20-200
GILSON	P-100 PIPETMAN	10-100
	P-200 PIPETMAN	50-200
HAMILTON	SOFTGRIP	10-100
	SOFTGRIP	30-300
LABNET	VE-100 LAB3 LABPETTE DIGITAL	10-100
	VE-200 LAB3 LABPETTE DIGITAL	20-200
NICHIRYO	200V-M TRIPLE VOLUME	20/100/200
OXFORD	BENCHMATE	40-200
RAININ	R-100 PIPET-PLUS	10-100
	R-200 PIPET-PLUS	50-200
SOCOREX	825 ACURA	10-100
	825 ACURA	20-200
	822 CALIBRA	10-100
	822 CALIBRA	20-200
TITERTEK	DIGITAL	40-200
VISTALAB	OVATION	10-100
	MLA	50-200