

## Adhesive Plate Seals

**Description:** A polyvinylidene chloride (PVdC) coated polypropylene seal with an acrylic adhesive. The film offers a convenient and cost effective adhesive cover for long and short-term storage. The non-adhesive tab allows quick removal and the film also creates a good reseal on top of a pierced foil seal.

**Ordering Information:** AB-0580 Adhesive Plate Seals 100 sheets

**Physical Properties:** Seal thickness: 55µm (without release liner)  
Temperature range: -20°C to 80°C

**Note:** This product is **not** recommended for high temperature applications such as PCR or for storing at temperatures below -20°C.

**For Research Purposes Only**