

2. SELECT YOUR DETECTION CHEMISTRY

1. SELECT YOUR CYCLER

QPCR Cyclers	P R O B E					S Y B R [®] G R E E N					
	Mix + ROX Vial	ROX Mix	Low ROX Mix	Mix (No ROX)	Capillary Mix (No ROX)	Mix + ROX Vial	ROX Mix	Low ROX Mix	Mix (No ROX)	Capillary Mix (No ROX)	Fluorescein Mix (No ROX)
ABI											
7000		●					●				
7300		●					●				
7500			●					●			
7700		●					●				
7900		●					●				
7900HT		●					●				
StepOne™		●					●				
StepOne™ Plus		●					●				
Bio-Rad											
iCycler™	●			●							●
MyiQ™	●			●							●
iQ™ 5	●			●		●			●		
Opticon™	●			●		●			●		
Opticon™ 2	●			●		●			●		
Chromo 4™	●			●		●			●		
MiniOpticon	●			●		●			●		
CFX96	●			●		●			●		
CFX384	●			●		●			●		
Cepheid											
SmartCycler™	●			●		●			●		
Corbett											
Rotor-Gene™	●			●		●			●		
Eppendorf											
Realplex	●			●		●			●		
Roche											
Lightcycler® 2.0					●					●	
Lightcycler® 480	●			●		●			●		
Stratagene											
Mx4000™	●		●			●		●			
Mx3000P™	●		●			●		●			
Mx3005P™	●		●			●		●			
Techne											
Quantica™	●			●		●			●		



3. CHOOSE YOUR MIX

QPCR	Absolute											
	200 x 25µl (2 x 1.25ml)	AB-1136/A	AB-1138/A	AB-1318/A	AB-1132/A	AB-1283/A	AB-1166/A	AB-1162/A	AB-1322/A	AB-1158/A	AB-1285/A	AB-1219/A
	1,600 x 25µl (16 x 1.25ml)	AB-1136/B	AB-1138/B	AB-1318/B	AB-1132/B	AB-1283/B	AB-1166/B	AB-1162/B	AB-1322/B	AB-1158/B	AB-1285/B	AB-1219/B
	400 x 25µl (1 x 5ml)	AB-1137/A	AB-1139/A	AB-1319/A	AB-1133/A	—	AB-1167/A	AB-1163/A	AB-1323/A	AB-1159/A	—	AB-1220/A
	4,000 x 25µl (10 x 5ml)	AB-1137/B	AB-1139/B	AB-1319/B	AB-1133/B	—	AB-1167/B	AB-1163/B	AB-1323/B	AB-1159/B	—	AB-1220/B
	Absolute Blue											
	200 x 25µl (2 x 1.25ml)	AB-4136/A	AB-4138/A	AB-4318/A	—	—	AB-4166/A	AB-4162/A	AB-4322/A	—	—	AB-4219/A
	1,600 x 25µl (16 x 1.25ml)	AB-4136/B	AB-4138/B	AB-4318/B	—	—	AB-4166/B	AB-4162/B	AB-4322/B	—	—	AB-4219/B
	400 x 25µl (1 x 5ml)	AB-4137/A	AB-4139/A	AB-4319/A	—	—	AB-4167/A	AB-4163/A	AB-4323/A	—	—	AB-4220/A
	4,000 x 25µl (10 x 5ml)	AB-4137/B	AB-4139/B	AB-4319/B	—	—	AB-4167/B	AB-4163/B	AB-4323/B	—	—	AB-4220/B
	Absolute Fast											
	200 x 25µl (2 x 1.25ml)	AB-4325/A	AB-4327/A	AB-4329/A	—	—	—	—	—	—	—	—
	1,600 x 25µl (16 x 1.25ml)	AB-4325/B	AB-4327/B	AB-4329/B	—	—	—	—	—	—	—	—
	400 x 25µl (1 x 5ml)	AB-4326/A	AB-4328/A	AB-4330/A	—	—	—	—	—	—	—	—
	4,000 x 25µl (10 x 5ml)	AB-4326/B	AB-4328/B	AB-4330/B	—	—	—	—	—	—	—	—
	Verso 1-Step											
200 x 25µl (2 x 1.25ml)	AB-4100/A	AB-4101/A	AB-4102/A	—	—	AB-4104/A	AB-4105/A	AB-4106/A	—	—	AB-4107/A	
400 x 25µl (1 x 5ml)	AB-4100/C	AB-4101/C	AB-4102/C	—	—	AB-4104/C	AB-4105/C	AB-4106/C	—	—	AB-4107/C	
Verso 2-Step												
200 x 25µl (2 x 1.25ml)	AB-4109/A	AB-4110/A	AB-4117/A	—	—	AB-4112/A	AB-4113/A	AB-4116/A	—	—	AB-4114/A	
400 x 25µl (1 x 5ml)	AB-4109/B	AB-4110/B	AB-4117/B	—	—	AB-4112/B	AB-4113/B	AB-4116/B	—	—	AB-4114/B	

Need help? Contact your local technical support team:

- United States +1 (800) 235-9880, option #4
- Rest of the world +44 (0)1372 723456, option #3

abgeneus.techservice@thermofisher.com
abgene.techsupport@thermofisher.com