

The Thermo Scientific Abgene SmartCap Mat rapidly seals an entire rack of Thermo Scientific storage tubes, ensuring reliable sealing from ambient temperatures down to $-80\text{ }^{\circ}\text{C}$.

Thermo Scientific Abgene® SmartCap Mat



Efficient Storage. The Thermo Scientific Abgene SmartCap Mat seals 96 samples at one time in a high-density storage format, affording maximum space saving benefits. Suitable for use with Thermo Scientific storage tubes ranging from 0.3 ml to 1.4 ml, the Abgene SmartCap Mat maintains a consistent seal down to $-80\text{ }^{\circ}\text{C}$.

Enhanced Materials. Abgene SmartCap Mats are made of an advanced thermoplastic elastomer that offers excellent chemical resistance and low vapor transmission rates, minimizing evaporation and water uptake in DMSO-based samples. To allow for sterilization of individual mats, the Abgene SmartCap Mat backing is autoclavable.

Ease of Installation and Access. Samples can be accessed in full rows of eight tubes or individually with the Thermo Scientific Abgene SmartCap Removal Tool, to effortlessly remove or re-install the caps. Sealing is secure and consistent every time when applied with the Thermo Scientific SuperSealer; after the SuperSealer applies the mat, simply pull up the backing sheet to reveal 96 individually sealed tubes.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables and can demonstrate how Abgene SmartCap Mats are an ideal solution for your storage tube sealing needs.

Thermo Scientific Abgene SmartCap Mat

Item No.	Description	Qty
AB-1413	Abgene SmartCap Mat	10 mats of 96 plugs/case
AB-1414	Abgene SmartCap Removal Tool	Each

Compatible Products

Thermo Scientific SuperSealer

Thermo Scientific Storage Tubes

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com
Europe: Tel: +44 (0) 161 486 2110 | email: matrix.eu.info@thermofisher.com
Asia: email: matrix.ap.info@thermofisher.com

www.thermo.com/matrix

Thermo
SCIENTIFIC