

dNTP Mix (High Concentration)

- Description:** dNTP mix (high concentration) containing four deoxynucleotides in a premixed solution, ready for immediate use.
- Ordering Information:** AB-1124 dNTP Mix (high concentration) – 1 x 1ml
25mM of each dNTP
- Concentration:** Combined 25mM dNTP mix containing 25mM of each dNTP (dATP, dCTP, dGTP, dTTP) at pH 7.0. The total dNTP concentration is 100mM.
- Storage Conditions:** Store at -20°C for up to 1 year. Shipped on ice within the UK and on dry ice for international and within the US.
- Application:** For primer extension procedures, add 0.4µl of dNTP mix in a total 50µl reaction mix to give a final concentration of 200µM of each dNTP.

For Research Purposes Only