

## Mineral Oil

**Description:** Mineral oil is a light, white oil suitable for overlaying molecular biology reactions such as PCR. It is supplied DNase, RNase and Protease free.

**Ordering Information:** AB-0262 Mineral Oil (1 vial of 30ml)

**Storage Conditions:** Stored and shipped at ambient temperature.

**Density:** 0.84g/ml

**Application:** Sample overlay.  
Not suitable for ABgene Diamond 384-well PCR Plate (AB-1111).

**For Research Purposes Only**