

Our Revco laboratory refrigerators and freezers are an economical choice for common non-critical laboratory requirements. Models provide dependable refrigeration systems and quality cabinet construction with basic control features.

Thermo Scientific Revco General Purpose Refrigerators, Freezers and Combinations



Refrigerators and Combinations

- Refrigerator operating temperature: +2°C to +10°C
- Freezer operating temperature: -12°C to -18°C
- High-strength shelves with interior door shelving for added inventory capacity
- Hermetically-sealed compressor for quiet operation, CFC-free refrigerants
- Built-in condensate drain on 20.7 cu. ft. models
- Magnetic door gaskets for positive seal; 20.7 cu. ft. models feature keyed door lock
- Power cord and plug provided for all models with 115V, 60Hz operation; NEMA 5-15R/5-15P
- Durable, white acrylic exterior finish and seamless interior

Laboratory Freezers

- Adjustable operating temperature: -12°C to -20°C
- High-strength shelves plus door shelving for increased storage capacity
- Hermetically-sealed compressor for quiet operation, CFC-free refrigerants
- High-density urethane foam insulation
- Magnetic door gaskets for positive seal; 20.7 cu. ft. models feature keyed door lock
- Power cord and plug provided for all models with 120V, 60Hz operation, NEMA 5-15R/5-15P
- Durable, white acrylic exterior finish and seamless interior

Thermo Scientific Revco General Purpose Refrigerators

Operating Temperature Range: +2°C to +10°C (manual defrost)

Model No.	Cu.Ft. (L)	Voltage (Hz)	Amps/ Breaker (Plug)	Shelves - Interior (Shelves - Door)	Interior Dimensions* H x D x W inches (mm)	Exterior Dimensions* H x D x W inches (mm)	Ship Weight lb (kg)
RLR2111A	20.7 (586)	120 (60)	5/15 (5-15R/5-15P)	4 (6)	57.5 x 23.5 x 27.0 (1460 x 600 x 690)	70.2 x 30.2 x 32 (1800 x 770 x 810)	250 (114)

Thermo Scientific Revco General Purpose Refrigerator/Freezer Combinations

Operating Temperature Range: +2°C to +10°C, Freezer -12°C to -18°C (manual defrost)

Model No.	Cu.Ft. (L)	Voltage (Hz)	Amps/ Breaker (Plug)	Shelves - Interior (Shelves - Door)	Interior Dimensions* H x D x W inches (mm)	Exterior Dimensions* H x D x W inches (mm)	Ship Weight lb (kg)
RCRF112A	8.6/244 (ref) 2.4/68 (frz)	120 (60)	3.8/15 (5-15R/5-15P)	4 (6)	ref: 35.3 x 19 x 20 (900 x 480 x 510) frz: 13 x 18 x 17.5 (330 x 460 x 440)	58.8 x 25.8 x 23.5 (1490 x 660 x 600)	182 (83)
RCRF252A	17.5/496 (ref) auto defrost 7.1/201 (frz)	120 (60)	5.7/15 (5-15R/5-15P)	4 (6) ref 1 (2) frz	ref: 37.8 x 24.5 x 32 (960 x 620 x 810) frz: 18.8 x 20.5 x 31 (480 x 520 x 790)	68.6 x 33.6 x 36 (1740 x 850 x 910)	142 (64)

Thermo Scientific Revco General Purpose Freezers

Operating Temperature Range: -12°C to -18°C (manual defrost)

Model No.	Cu.Ft. (L)	Voltage (Hz)	Amps/ Breaker (Plug)	Shelves - Interior (Shelves - Door)	Interior Dimensions* H x D x W inches (mm)	Exterior Dimensions* H x D x W inches (mm)	Ship Weight lb (kg)
REF417A	5.0 (141)	115 (60)	5/15 (5-15R/5-15P)	2 (3)	29 x 19 x 20 (737 x 476 x 495)	34 x 21 x 27 (851 x 540 x 673)	84 (38)
REF1617A	17.3 (490)	120 (60)	1.8/15 (5-15R/5-15P)	3 (5)	51.5 x 21.5 x 27 (1310 x 550 x 690)	65.2 x 28 x 32 (1660 x 710 x 810)	235 (107)
REF2117A	20.7 (586)	120 (60)	1.8/15 (5-15R/5-15P)	4 (6)	57 x 23 x 27 (1450 x 580 x 690)	70.2 x 30.2 x 30 (1780 x 770 x 810)	250 (114)

*All dimensions and specifications subject to change without notice. External dimensions include allowances for cabinet door, hinges, latches, handles, and control panel protrusions. †No plug provided; requires hardwire connection by customer.



© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

www.thermo.com/cold

North America: USA/Canada +1 866 984 3766

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2000, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254, Netherlands +31 76 571 4440, Nordic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India +91 22 6716 2200, Japan +81 45 453 9220, Other Asian countries +852 2885 4613

Countries not listed: +49 6184 90 6940 or +33 2 2803 2000

DS-LECS-REVCOREF-0808

Thermo
SCIENTIFIC