

Blood Gas Supplies for use on IL 1610 and 1620

Products:

Calibration Buffers
 Flush Solution
 Calibration Gases
 Thermal Printer Papers

Shelf Life:

2 year outdate
 2 year outdate
 3 year outdate
 5 year outdate

Intended Use:

6.840/7.384 Buffers To provide calibration points for the pH electrode on the above instruments.
 Flush Solution Kit To flush the sample flow path for the above instruments.
 Calibration Gases To provide calibration points for the pCO₂ and pO₂ electrodes on the above instruments.
 Printer Papers To print test results.



For Use On:

Instrumentation Laboratory 1610
 Instrumentation Laboratory 1620

Analytes:

pH, pCO₂, pO₂
 pH, pCO₂, pO₂

<i>OEM #</i>	<i>Thermo ALKO #</i>	<i>Product Description</i>	<i>OEM Pkg</i>	<i>Thermo ALKO Pkg</i>
31060	A310-60	6.840 pH Calibration Buffer	4x500 ml	4x500 ml
33360	A333-60	7.384 pH Calibration Buffer	4x500 ml	4x500 ml
33352	A333-52Z	Flush Solution Kit	4x500 ml	4x500 ml
74389	A743-89	Printer Paper, 1610	6/pkg	6/pkg
190016	A800-75	Printer Paper, 1620/30/40/50	4/pkg	4/pkg
52321	A523-21	Gas 5% CO₂ 20% O₂ Bal N₂	20 scf	20 scf
52220	A522-20	Gas 10% CO₂ Bal N₂	20 scf	20 scf

Additional Information:

1. Blood Gas/Electrolyte Controls (A701 series) have these assay values.
2. Combination Quality Controls (A100 series) are also available with assays for these instruments.