

# Thermo ALKO

## Blood Gas Supplies for use on AVL 940, 945, 990, 995 and 995 Hb

Correlation data is available upon request.

### Shelf life

Most Reagents: 4 yrs from date of manufacture

Calibration Gases: 3 yrs

Quality Controls: 4 yrs from date of manufacture



<b>OEM #</b>	<b>Thermo ALKO #</b>	<b>Product Description</b>	<b>OEM Pkg</b>	<b>Thermo ALKO Pkg</b>
BP0136	A007-7	7.383 Type 1 Buffer	3 x 90 mL	3 x 90 mL
BP0137	A007-8	6.841 Type 2 Buffer	3 x 90 mL	3 x 90 mL
BP9076	A907-6	Rinse Additive Solution	100 x 1 mL	250 mL
BP1025	A009-9A	Cleaning Solution A	100 mL	250 mL
BP1026	A102-6	Cleaning Solution B	100 mL	250 mL
HC0014	A011-0	pO <sub>2</sub> Electrolyte Solution	100 mL	250 mL
HC0010	A011-1	pCO <sub>2</sub> Electrolyte Solution	100 mL	250 mL
BP0521	A477-009Z	Deproteinizing Cleaning Sol	100 mL	6 x 125 mL
BP0134	A012-1	pH Reference Solution	3 x 90 mL	3 x 90 mL
n.a.	A300-308-A	Gas 99.5% CO <sub>2</sub>	n.a.	48 scf
HC9027	A701-001	ALKontrol, Level 1	30 x 1.7 mL	30 x 2.5 mL
HC9028	A701-002	ALKontrol, Level 2	30 x 1.7 mL	30 x 2.5 mL
HC9029	A701-003	ALKontrol, Level 3	30 x 1.7 mL	30 x 2.5 mL
n.a.	A701-123	ALKontrol, Tri-Level	n.a.	30 x 2.5 mL
HC9030	A500-004	ALKontrol, Level 4 High O <sub>2</sub>	30 x 1.7 mL	30 x 2 mL

### Products:

7.383 Type 1 Buffer

6.841 Type 2 Buffer

Rinse Additive Solution

Cleaning Solution A

pO<sub>2</sub> Electrolyte Sol

pCO<sub>2</sub> Electrolyte Sol

Deproteinizing Sol

pH Reference Sol

### Intended Use:

A calibration point for pH electrode of AVL pH/Blood Gas analyzers.

A calibration point for pH electrode of AVL pH/Blood Gas analyzers.

To rinse the sample flow path of AVL pH/Blood Gas analyzers.

To clean the sample flow path of pH/Blood Gas and Electrolyte analyzers.

To fill the pO<sub>2</sub> electrode of AVL pH/Blood Gas analyzers.

To fill pCO<sub>2</sub> electrode of AVL pH/Blood Gas analyzers.

To clean the sample flow path of pH/Blood Gas and Electrolyte analyzers.

Provide constant electrical potential through the reference electrode of AVL pH/Blood Gas analyzers.

The names of various manufacturers, their instruments, and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Thermo ALKO expressly disclaims any affiliation with them or sponsorship by them.