

Thermo ALKO

Thermo ALKO ISE Supplies for use on Beckman Lablyte 800, 810, 820, 830



OEM #	Thermo ALKO #	Product Description	Instrument Models	OEM Pkg	Thermo ALKO Pkg
443930	A093-6	Standard A	800,810,830	3 x 430 mL	3 x 430 mL
443935	A084-0	Standard A	820	3 x 430 mL	3 x 430 mL
443940	A800-094A	Standard B	800,810	3 x 90 mL	3 x 90 mL
443945	A084-1	Standard B	820	3 x 90 mL	3 x 90 mL
443950	A093-7	Standard B	830	3 x 90 mL	3 x 90 mL
443955	A093-8	Standard C	830	3 x 90 mL	3 x 90 mL
443965	A800-095	Reference Solution	800,810,820,830	3 x 90 mL	3 x 90 mL
n.a.	A381-013	Na/K/Li Electrode Fill Sol	820,830	n.a.	2 x 125 mL
n.a.	A214-0	Na/K/Cl Electrode Fill Sol	800,810	n.a.	2 x 125 mL
443970	A000-18	Na/pH Electrode Cond Sol	800,810,820,830	100 mL	250 mL
443980	A009-9A	Cleaning Solution A	800,810,820,830	100 mL	250 mL
756461	A041-3	Na Electrode	800,810,820,830	each	each
756462	A035-9	K Electrode	800,810,820,830	each	each
756467	A071-3	Reference Electrode	800,810,820,830	each	each
756466	A084-4	Reference Housing	800,810,820,830	each	each
949223	A820-091	Thermal Printer Paper	800,810,820,830	5 rolls	5 rolls
946420	A821-111A	Pump Tubing	800,810,820,830	1 set	2 sets
756468	A105-400	Intake Tubing	800,810,820,830	1 set	1 set
n.a.	A701-001	ALKontrol, Level 1	830	n.a.	30 x 2.5 mL
n.a.	A701-002	ALKontrol, Level 2	830	n.a.	30 x 2.5 mL
n.a.	A701-003	ALKontrol, Level 3	830	n.a.	30 x 2.5 mL
443985	A701-123	ALKontrol, Tri-Level	830	12 x 1 mL	30 x 2.5 mL

Additional Information:

1. Instrument parameters

800	Na/K	820	Na/K/Ca/pH
810	Na/K/Cl	830	Na/K/Li

2. Shelf Life

Most Reagents:	4 yrs from date of manufacture
Electrodes and Accessories:	Inquire for specific outdates
Quality Controls:	4 yrs from date of manufacture

The names of various manufacturers, their instruments, and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Thermo ALKO expressly disclaims any affiliation with them or sponsorship by them.