

GRAMS 8.0 Quick Start Installation Instructions

Please note that you must have Administrative rights on your computer in order to install successfully.

1. Visit <http://gramssuite.com/downloads/grams8> to reach the GRAMS download page.
2. After performing the first two steps on the download page, continue with Step #3 here (Step #3 on the site is simply a pointer to these instructions.)
3. Open Windows Explorer.
4. In Windows Explorer, browse to the location of the installation file (GRAMS80Install.exe). Double-click this file to start the installation.
5. Click "Unzip" in the WinZip Self-Extractor dialog. The files will be unzipped to the folder C:\GRAMSCDImage.
6. When unzip is finished, click the "OK" button. The install will be started automatically.
7. Click "Install Products" on the "GRAMS Suite Installer" screen.
8. Read the license agreement. Press the "Accept" button if you accept the terms of the license agreement.
9. On the "Open License File" screen, press the "Browse" button to browse to the location of your license.txt file for GRAMS 8. (If you purchased a GRAMS product, this file was Emailed to you. If you are trying the demonstration copy, this is the file you downloaded in Step #1 from the download page.) Select the license.txt file and press "Open".
10. Press the "Next" button to continue the installation.
11. The "Installation Folder" should have a default path of "C:\Program Files\Galactic". We recommend that you do not change this path.
12. Press the "Typical" button to continue the installation. The installation should proceed automatically.
13. If the installation was successful, you will get a message that says "Finished installing X products" (where X is the number of products you have purchased). Press the OK button.
14. Press the "Exit" button to exit the installation program.
15. Before running the GRAMS suite, you must install the "client" portion of the software for each application. Click on your Windows Start menu and click "Programs" and then select "Thermo GRAMS Suite" and then click on the short cut to one of the GRAMS Suite Applications (GRAMS AI, GRAMS3D, Spectral ID, and Spectral DB). You will need to repeat this for each of the application listed above that you have purchased.
16. (Optional) Once you are sure the installation has been successful, you may remove the temporary folder: C:\GRAMSCDImage.