

Get Mercury Out of Your Laboratory

If you're working with TRIS buffers or samples containing proteins and enzymes, chances are you're using a Calomel pH electrode which contains dangerous mercury!

Non-mercury

Orion Mercury Free Electrodes



- Does your lab have restrictions on mercury-containing products?
- Are you unknowingly introducing mercury into your lab through your pH electrodes?
- Are you concerned with mercury clean-up if your pH electrode breaks?

Why handle hazardous material when you have a choice.

Thermo offers non-mercury Orion electrodes specifically developed to outperform even in your most difficult samples. Choose from a variety of options from a value positioned Ag/AgCl electrode or a AquaPro where low maintenance is important or select the best of the best- ROSS Ultra® refillable with unmatched reliability, performance, plus a 2 year warranty.

ROSS and the COIL traddress are trademarks of Thermo Electron Corporation. US Patent 6,793,787

Thermo offers non-mercury Orion electrodes specifically developed to outperform even in your most difficult samples. Plus, Orion electrodes with APWP, BNUWP or DJWP connections will work on any pH meter with a BNC connector.

GOOD

Value positioned double junction Ag/AgCl pH electrodes with 0.02 to 0.05 pH unit precision and a 1 year warranty

9110DJWP	Semi-micro glass body, refillable pH electrode, volumes down to 0.2mL
9102DJWP	Standard glass body, refillable pH electrode
9156DJWP	Standard epoxy body, non-refillable pH electrode*

BETTER

AquaPro non-refillable polymer pH electrodes are low maintenance yet provide fast response times with 0.02 pH unit precision and a 1 year warranty

9103APWP	Semi-micro glass body pH electrode, volumes down to 0.2mL
9104APWP	Rugged bulb standard glass body pH electrode
9115APWP	Semi-micro epoxy body pH electrode, volumes down to 0.2mL
9156APWP	Standard epoxy body pH electrode

BEST

ROSS Ultra® refillable pH electrodes provide unmatched stability and reliability with the fastest response, highest 0.01 pH unit precision and a 2 year warranty.

8102BNUWP	Standard glass body pH electrode
8156BNUWP	Standard epoxy body pH electrode
8103BNUWP	Semi-micro glass body pH electrode, volumes down to 0.2mL
8115BNUWP	Semi-micro epoxy body pH electrode, volumes down to 0.2mL

* 9156DJWP Gel filled electrode is warranted to be free from defects in material and workmanship for a period of six (6) months from date of purchase by customer or twelve (12) months from date of shipment from Thermo Electron, whichever is earlier.

ROSS and the COIL tradenames are trademarks of Thermo Electron Corporation. US Patent 6,793,787