

<b>PART #</b>	<b>DESCRIPTION</b>
7410	BENCH, OPTICAL
9931	BUTTON, ARROW SIDE
9932	BUTTON, ARROW UP/DOWN
9928	BUTTON, ENTER
9930	BUTTON, HELP
9929	BUTTON, MENU
9927	BUTTON, RUN
9958	CABLE ASSY. ( HC 11 )
6279	CABLE RS232 ASSY 6'
7356	CAPACITOR (4762 MOTOR) 220 VOLTS
7336	CAPILLARY, 18 MIL L. (BLUE/PURPLE)
8938	DETECTOR PRE-AMP (SPECIFY 48C/48CHL)
8630	FAN GUARD ASSEMBLY (W/Foam) "C-SERIES"
4299	FAN, COOLING, SPRITE 115V/50/60
4304	FEET, RUBBER (PK/4) ("C" Series Instrument Case)
8544	FEET, RUBBER (SHOCK MOUNTS) (Optical Benches)
7380	FILTER ASSEMBLY, BANDPASS
7358	FILTER WHEEL, CORRELATION (CO)
9934	FLOW TRANSDUCER ASSY. ( HC11 )
4508	FUSE 1 1/2 AMP S/B (PK/5)
4510	FUSE 3 AMP S/B (PK/5)
5924	HEATER RESISTOR, 50 OHM (115 V) (Specify Instrument)
9886	INSTRUMENT, COVER
9998	MANUAL, MODEL 48C
7413	MIRROR, ENTRANCE (48)
7414	MIRROR, EXIT (48)
7412	MIRROR, FIELD (48) (6 Holes)
7411	MIRROR, RELAY (48) (3 Holes)
4735	MOTOR BODINE ( 48C )
7211	OPTICAL BENCH, COMPLETE WITH FILTER WHEEL (48C)
4807	O-RING (FOR RELAY/FIELD MIRRORS)
4827	O-RING (FOR THE FILTER INTERFERENCE ASSEMBLY)
4818	O-RING (PK/10) (FOR ENTRANCE MIRROR )
4826	O-RING (PK/5) (FOR EXIT MIRROR)
8331	PCB A/D (48C)
9843	PCB C-LINK (Must Specify Program Version)
9835	PCB CLOCK
9904	PCB CONNECTOR INTERFACE
9839	PCB DIGITAL/ANALOG
11220	PCB DISPLAY
9989	PCB INPUT ( 48C )

<b>PART #</b>	<b>DESCRIPTION</b>
9831	PCB MOTHER ( 48C )
8931	PCB POWER SUPPLY ( 48C )
8933	PCB POWER SUPPLY ( BIAS ) ( 48C )
9837	PCB PROCESSOR (Must Specify Program Version)
9950	PCB PUSH BUTTON ( B1 )
9952	PCB PUSH BUTTON ( B2 )
8935	PCB TEMP. CONTROLLER ( 48C )
11519	PLUG 8 PIN (Rear Panel Terminal Strip)
14010	PLUG RECEPTACLE (49C)
8165	POWER CORD, 115V
8709	PUMP (115V/50/60Hz) Specify Voltage and Model Number)
8550	PUMP KNF 115V/60HZ (PU603N056)
8551	PUMP KNF 230V/50HZ
8606	PUMP REPAIR KIT, KNF (FOR P/N 8550)
11289	RECEPTACLE 2 POLE
18061	RUBBER FEET (Insulator) (Feet for the Pumps)
7360	SENSOR ASSEMBLY, OPTICAL
7374	SENSOR OPTICAL
7384	SOURCE COVER
7361	SOURCE IR
7385	SOURCE MOUNT IR
12118	SWITCH POWER MODE SP3T
5382	THERMISTOR IM-1002-A5
9877	TRANSDUCER, PRESSURE ASSY.
5496	TUBING, SILASTIC (FOOT)
7368	VALVE SOLENOID 3-WAY