

PART #	DESCRIPTION
11526	BATTERY CELL (Panasonic)
81P751	BIAS POGO ASSY
81P737	BYPASS BULKHEAD ASSY
6279	CABLE RS232 ASSY 6'
81P717	COLLECTOR ASSY
16004	CONTROLLER AIRFLOW (Make-up Air Solenoid)
18186	DIFFUSER
81P719	FID SAMPLE FRIT ASSY
12190	FILTER 0.5V
4070	FILTER ELEMENT, FAN
16011	FLAT WASHER
14007	FUSE 1.5 AMP 250V AC SB (PK/5)
14011	FUSE, TRANSFORMER
10260	GAUGE PANEL MOUNT
15018	IGNITOR ASSEMBLY
81P534	INSULATOR (BASE)
81P537	INSULATOR (UPPER DETECTOR)
81P716	JET ASSY
11032	O-RING (PK/2)
6646	O-RING VITON (PK/5)
16107	PCB ASSY, A/D
16077	PCB ASSY, ANALOG BOARD
16106	PCB ASSY, MOTHER BOARD
16072	PCB ASSY, POWER SUPPLY
16083	PCB ASSY, REAR CONNECTOR
9843	PCB C-LINK (Must Specify Program Version)
11220	PCB DISPLAY
3513	PCB I/O CONTROLLER
9837	PCB PROCESSOR (Must Specify Program Version)
81P728	PRESSURE SENSOR ASSEMBLY
81P722	PUMP ASSY - Specify Voltage 110V or 220V
81P762	PUMP ASSY (51CHT)
81P759	PUMP CABLE ASSEMBLY
18070	PUMP REPAIR KIT, (FOR P/N 17012)
8606	PUMP REPAIR KIT, KNF (FOR P/N 8550)
14954	Regulator (Air or Fuel)
8658	REGULATOR (NOROGEN) Make-up Air Regulator
81P739	SIGNAL POGO ASSY
14951	SOLENOID VALVE 2-WAY (55C)
12118	SWITCH POWER MODE SP3T
81P729	THERMOCOUPLE ASSEMBLY

Thermo Electron Corporation
27 Forge Parkway
Franklin, MA 02038
Telephone # 866-282-0430
Fax # 508-520-2800

Thermo Electron
51C

PART #	DESCRIPTION
81P758	ZERO SPAN FRIT ASSEMBLY
16021	ZERO/SPAN SOLEDOID