

# Material Safety Data Sheet

**ThermoFisher**  
SCIENTIFIC

## NERL Diagnostics, LLC

Subsidiary of Thermo Fisher Scientific Inc.

14 Almeida Avenue

East Providence, RI 02914

Telephone: 1-800-556-7575

Fax: 1-401-438-2454

### I. PRODUCT INFORMATION

#### Negative Fecal Hemoglobin Control

	<u>Chemical Formula</u>	<u>CAS#</u>	<u>Threshold Limit Value</u>	
Hazardous Components:	Microcide III (blend)	n/a	26172-55-4/2682-20-4	n/a

### II. TOXICITY HAZARDS

Microcide III – Irritant to eyes and skin. May cause burns upon prolonged contact (corrosive).

### III. PHYSICAL DATA

<b>Appearance:</b>	Brown liquid	<b>Evaporation Rate:</b>	N/A
<b>Odor:</b>	N/A	<b>Specific Gravity:</b>	N/A
<b>Boiling Point:</b>	N/A	<b>pH:</b>	N/A
<b>Osmolality:</b>	N/A	<b>Water Solubility:</b>	Soluble

### IV. FIRE & EXPLOSION INFORMATION

<b>Flash Point:</b>	N/A
<b>Fire Extinguishing Media:</b>	Dry chemical, carbon dioxide, water spray, or foam.
<b>Fire Fighting Procedures:</b>	Cool fire-exposed areas with water spray.
<b>Explosion Hazards:</b>	N/A

### V. HEALTH HAZARDS

<b>Threshold Limit Value:</b>	Unknown
<b>Effects of Overexposure:</b>	N/A
<b>Routes of Exposure:</b>	Ingestion: May be harmful if swallowed. May be toxic if ingested.
	Inhalation: N/A

