

Fluoride measurements from On-Line to the Lab have never been easier... introducing the NEW Orion 2109XP Fluoride Monitor.

For a limited time purchase an Orion 2109XP monitor and receive a FREE Orion 4-Star Fluoride Benchtop Meter Kit! \*

**\$ 2,148.00 Savings!**

**Thermo Scientific  
Orion 2109XP  
Fluoride Monitor**



**Orion 2109XP Performance Benefits**

The Orion 2109XP fluoride monitor offers accuracy, reliability and ease-of-use to best meet the demands for fluoride analysis with complete assurance. Validation of measurements is fast and simple using the Orion 2109XP and the Orion 4-Star fluoride benchtop ISE system. Create confidence in your water quality results when it matters the most!

- Detection limit of 10 ppb - continuous fluoride monitoring in a wide variety of samples using the Orion 2109XP
- Based upon EPA method for fluoride using an ion selective electrode (ISE)
- Fastest, most stable measurements - reduce calibration cycles due to drift with superior Orion fluoride electrodes

- Eliminate problems with turbidity and suspended solids in samples - the Orion 2109XP provides true accuracy in measuring total fluoride
- Easy to use while maximizing uptime - simple step-by-step scrolling instruction for setup, calibration, operation and diagnostics menus
- Measurements at a glance from any distance using the large operator friendly visible backlit display
- Minimize operator time and maintenance without use of complicated moving parts or pumps that often require frequent attention and expensive spare part expenditures

\*Limited time only, while supplies last, offer commencing March 10, 2008 and is valid for US customers only. This offer may not be combined with any other coupon, discount, or package. This promotion must be mentioned to receive the offer. This offer applies only for the products listed and purchased after March 10, 2008. All prices are in U.S. dollars. This offer may be discontinued at any time and is void where prohibited by law or company policy. Offer expires December 31, 2008. Please allow 4-6 weeks for delivery of promotional product.

**Promotion valid for US direct customers only.**

## Ordering Information

Simply place your order for the Orion 2109XP fluoride monitor package of your choice and Orion Cat. No. 1010103, an Orion 4-Star fluoride benchtop meter kit will be shipped with your 2109XP monitor at no charge.

- Promotion begins March 10, 2008 and is valid through December 31, 2008

Cat. No.	Fluoride Monitor Packages
<b>2109XPEN</b>	<b>Fluoride Monitor with Protective Enclosure Package</b> – Includes fluoride monitor with protective enclosure, fluoride sensing electrode (100029), reference electrode (100056), ATC (2100TP), standards solution kit (180940), CD user guide and options kit.
<b>2109XP</b>	<b>Fluoride Monitor Instrument Only Package</b> – Includes fluoride monitor with fluoride sensing electrode (100029), reference electrode (100056), ATC (2100TP), standards solution kit (180940), CD user guide and options kit.
Cat. No.	FREE Meter Kit Offer
<b>1010103</b>	<b>Orion 4-Star Fluoride Benchtop Meter Kit</b> <ul style="list-style-type: none"> <li>• Orion 4-Star benchtop meter (for ISE/pH measurements)</li> <li>• 9609BNWP fluoride ion selective electrode</li> <li>• Stirrer probe</li> <li>• Fluoride standards and solutions (pint bottles)</li> <li>• Swing arm electrode stand</li> <li>• Fluoride in drinking water method note</li> </ul>

For more details regarding this promotion contact us at 978-232-6000 or 1-800-225-1480, or visit us on the web at [www.thermo.com/process](http://www.thermo.com/process).

