

ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : CALCIUM

Catalogue no. : 981367

Lot no. : E117

Package : 6 x 60 ml

Storage Conditions : 2...25 °C

Expiry date : 2011-01

Quality Control

Performance test:

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol /l)</u>
Abtrol (D390A)	3,13	2,94-3,18
Nortrol (D389A)	2,34	2,24-2,48

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 17 / 3 2009



Hannele Laine
Quality Assurance Specialist
Quality Regulatory and Compliance
Clinical Diagnostics Finland

Konelab 30i (QC)
 SN: P1118753
 Thermo Fisher Scientific (CDX Finland)

12.3.2009 10:29

981367 E117 2011-01


Test CA

Accepted -

Factor 6.385
 Bias 0.002

Coeff. of det. 0.999960

Errors



Calibrator	Response	Calc. con.	Conc.	Errors
1 Water	0.002	-0.000	0.000	
2 Water	0.002	0.000	0.000	
3 sCal	0.387	2.459	2.470	
4 sCal	0.390	2.481	2.470	
5 Abtrol(control)	0.485	3.088	3.060	
6 Nortrol(control)	0.374	2.378	2.360	