

ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : Magnesium

Catalogue no. : 981385

Lot no. : E321

Package : 12 x 20 ml

Storage Conditions : 2...25 °C

Expiry date : 2011-05

**Quality Control
Performance test:**

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol/l)</u>
Abtrol (D390)	1,35	1,18 - 1,38
Nortrol (D389)	0,80	0,73 - 0,85

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 13 / 7 2009

Marja Keurulainen
Marja Keurulainen
Quality Chemist
Quality Regulatory and Compliance
Clinical Diagnostics Finland

Konelab 60i (QC)
SN: A5120360

10.7.2009 14:30

081385 lot E321, 2011-05


Test MG

Accepted 10.7.2009 13:08

Factor 4.703
Bias -0.005

Coeff. of det. 0.999845

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	Water	-0.004	0.001	0.000	
2	Water	-0.005	-0.001	0.000	
3	sCal	0.212	1.021	1.030	
4	sCal	0.216	1.039	1.030	
5	Nortrol(control)	0.159	0.771	0.790	
6	Abtrol(control)	0.280	1.339	1.280	