

## ANALYSIS CERTIFICATE

Manufacturer Name and Address: Thermo Fisher Scientific Oy  
Ratastie 2, FIN-01620 Vantaa, Finland

Product Name : Magnesium

Catalogue no. : 981385

Lot no. : E478

Package : 12 x 20 ml

Storage Conditions : 2...25 °C


Expiry date : 2011-05

**Quality Control  
Performance test:**

<u>Control (lot)</u>	<u>Result (mmol/l)</u>	<u>Target range (mmol/l)</u>
Abtrol (D390)	1,33	1,22 - 1,44
Nortrol (D389)	0,82	0,75 - 0,89

This product with the mentioned lot number was analysed by our authorised quality control department and confirmed that it meets its specifications.

Vantaa 10 / 12 2009

  
Hannele. Laine  
Quality Assurance Specialist  
Quality Regulatory and Compliance  
Clinical Diagnostics Finland

Konelab 60i (QC)  
SN: A5120360

9.10.2009 10:58

981385

E478

2011-05


Test MG

Accepted 9.10.2009 10:25

Factor 4.803  
Bias -0.005

Coeff. of det. 0.999949

Errors



Calibrator	Response	Calc. con.	Conc.	Errors
1 Water	-0.004	0.0032	0.0000	
2 Water	-0.006	-0.0032	0.0000	
3 sCal	0.209	1.0259	1.0300	
4 sCal	0.210	1.0341	1.0300	
5 Nortrol (control)	0.158	0.7842	0.7900	
6 Abtrol (control)	0.271	1.3248	1.2800	