

ROTORS

TUBES, BOTTLES

(Number of Tubes per Rotor)

NOMINAL CAPACITY (ml)	MATERIAL ¹	STYLE ²	DIMENSIONS (O.D. x L mm)	CATALOG NUMBER ³	F-28/50	F-28/36	F-28/13	F-16/250	S-20/36	S-20/20	S-20/17
1	CAB	T	7 x 50	03103	40						
1.5	MICROTUBE	—	—	—	20	12	16				
2	CAB	T	7 x 75	03235	40						
3	PY	T	10 x 75	03100	20			30			
3.5	PA	UC	11 x 60	03905		12	16				
4	PC	FL	11 x 75	03104	20			72			
4	PP	FL	11 x 75	03105	20			72			
4	SG	—	10 x 75	—				72			
4.4	PC	TH	11 x 60	54361		12	16				
4.4	PA	TW	11 x 60	03955		12	16				
5	PY	T	12 x 75	03102				72			
5	SG	—	12 x 75	—				72			
6	PA	UC	13 x 66	03945		12	16				
6	PET	CC	13 x 66	01086		12	16				
7	PP	FL	13 x 100	03121	10						
7	PC	FL	13 x 100	03120	10						
10	PC	OR	16 x 83	03020	10	12	16				
10	PP	OR	16 x 83	03929	10	12	16				
10	PC	T	16 x 71	03237		12	16				
10	PA	T	16 x 71	03238		12	16				
10	PA	TW	13 x 100	03124	10						
10	PY	—	16 x 100	—	10						
11.5	PA	UC	16 x 76	03987		12	16				
11.5	PET	CC	16 x 76	01087		12	16				
12	PC	FL	16 x 100	03115	10			42			
12	PP	FL	16 x 100	03116	10			42			
12	PY	T	18 x 102	03109	10						
12.5	PA	TW	16 x 76	03128		12	16				
14	PY	T	18 x 120	03111				30			
15	CO	T	18 x 120	03286	10						
15	—	TC,CN	17 x 120	—				30			
15.5	PA	UC	16 x 100	03409							6
16	PC	FL	18 x 100	03243	10						
16	PP	FL	18 x 100	03244	10						
17	PA	TW	16 x 102	03126							6
20	PA	TW	25 x 51	03288						6	
30	PC	OR	25 x 94	03719	10	12					
30	PP	OR	25 x 94	03718	10	12					
30	CO	T	24 x 105	03163	10			18			

(1) MATERIAL CODES

- CAB—Cellulose Acetate Butyrate
- PY—Pyrex® Glass
- PA—Polyallomer
- PS—Polysulfone
- CO—Corex® Glass
- PC—Polycarbonate
- PP—Polypropylene
- SS—Stainless Steel
- PET—Polyethylene-terephthalate
- SG—Standard Glass

(2) STYLE CODES

- T—Thick-Walled
- UC—ULTRACRIMP®
- TH—Threaded
- DS—DRY-SPIN Bottle
- TW—Thin-Walled
- OR—Oak Ridge Bottle
- FL—Flanged
- CN—Conical
- G—Graduated
- TC—Tissue Culture
- CC—CLEARCRIMP®

⁽³⁾If catalog number is not supplied, tube/bottle is not available from SORVALL.

All trademarks are the property of their respective owners.

TO ORDER, CALL 1-800-522-SPIN (7746)

TUBES, BOTTLES **ROTORS**
(Number of Tubes per Rotor)

NOMINAL CAPACITY (ml)	MATERIAL ¹	STYLE ²	DIMENSIONS (O.D. x L mm)	CATALOG NUMBER ³	ROTORS						
					F-28/50	F-28/36	F-28/13	F-16/250	S-20/36	S-20/20	S-20/17
34	PA	UC	25 x 80	03469					6		
35	PA	UC	25 x 97	03989		12					
35	PA	CN, OR	29 x 105	—	10			6			
35	PC	CN, OR	29 x 105	—	10			6			
36	PA	TW	25 x 89	03141		12			6		
50	SS	T	28 x 101	00579	10						
50	SS	FL	28 x 101	00517	10						
50	PC	OR	28 x 103	03480	10						
50	PP	OR	28 x 103	03530	10						
50	PA	TW	28 x 104	03139	10						
50	PC	FL	29 x 102	03146	10						
50	PP	FL	29 x 102	03147	10						
50	PA	UC	29 x 108	03528	10						
50	—	TC, CN	30 x 115	—				6			
100	PY	B	44 x 137	03230				6			
125	CO	B	53 x 118	03231				6			
150	PC	T	45 x 132	03143				6			
150	CO	B	53 x 132	03229				6			
250	PP	DS	61 x 146	03069				6			
250	PC	DS	61 x 146	03098				6			
250	PC	OR	61 x 124	03939				6			
250	PP	OR	61 x 124	03937				6			
250	PA	OR	60 x 120	54477				6			

(1) MATERIAL CODES

- CAB—Cellulose Acetate Butyrate
- PY—Pyrex® Glass
- PA—Polyallomer
- PS—Polysulfone
- CO—Corex® Glass
- PC—Polycarbonate
- PP—Polypropylene
- SS—Stainless Steel
- PET—Polyethyleneterephthalate

(2) STYLE CODES

- T—Thick-Walled
- UC—ULTRACRIMP®
- TH—Threaded
- DS—DRY-SPIN Bottle
- TW—Thin-Walled
- OR—Oak Ridge Bottle
- FL—Flanged
- CN—Conical
- G—Graduated
- TC—Tissue Culture
- CC—CLEARCRIMP®
- B—Bottle

⁽³⁾If catalog number is not supplied, tube/bottle is not available from SORVALL.

All trademarks are the property of their respective owners.

TO ORDER, CALL 1-800-522-SPIN (7746)

B-1