

TUBES, BOTTLES FOR F-28/36 (continued)

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
10	6.8	Polyallomer Thick-Walled Tube	28,000	16 × 71	03238 (50)	—	—	02210 (4)	1	Adapter also available in pkg. of 1 00613 .
10	6.8	Polycarbonate Thick-Walled Tube	28,000	16 × 71	03237 (50)	—	—	02210 (4)	1	Adapter also available in pkg. of 1 00613 .
10	9.1	Polypropylene Oak Ridge Bottle	28,000	16 × 80	03929 (50)	03613 (25)	Nylon, Screw Cap	02210 (4)	1	Adapter also available in pkg. of 1 00613 . Maximum speed is 21,000 rpm when partially filled.
10	9.7	Polycarbonate Oak Ridge Bottle	28,000	16 × 80	03020 (25)	03613 (25)	Nylon, Screw Cap	02210 (4)	1	Adapter also available in pkg. of 1 00613 .
6	5.7	Polyallomer ULTRACRIMP® Tube	28,000	13 × 66	03945 (50)	Included with tubes	Plug and Aluminum Cap	02209 (4)	1	Requires Rotor Cap 03992 , Tube Racks 03921 , and ULTRACRIMP® and CLEARCRIMP® Sealing Tool 03920 . Adapter also available in pkg. of 1 00612 .
6	5.7	Polyethylene-terephthalate CLEARCRIMP® Tube	28,000	13 × 66	01086 (50)	Included with tubes	Plug and Aluminum Cap	02209 (4)	1	Requires Rotor Cap 03992 , Tube Racks 03921 , and ULTRACRIMP® and CLEARCRIMP® Sealing Tool 03920 . Adapter also available in pkg. of 1 00612 .
4.4	3.3	Polyallomer Thin-Walled Tube	28,000	11 × 60	03955 (50)	54366 (1)	Multipiece, Selfseal	02208 (4)	1	Requires Tube Cap Vise 52176 and Tube Removal Tool 11538 . Adapter also available in pkg. of 1 00611 .
4.4	2.7	Polycarbonate Threaded Thick-Walled Tube	28,000	11 × 60	54361 (50)	54365 (1)	Aluminum, Screw Cap	02208 (4)	1	Requires Tube Cap Vise 54377 and Tube Removal Tool 11538 . Adapter also available in pkg. of 1 00611 .
3.5	3.5	Polyallomer ULTRACRIMP® Tube	28,000	11 × 60	03905 (50)	Included with tubes	Plug and Aluminum Cap	02208 (4)	1	Requires Tube Cap 03091 , Tube Racks 03773 , and ULTRACRIMP® and CLEARCRIMP® Sealing Tool 03920 . Adapter also available in pkg. of 1 00611 .
1.5	—	Microtube	—	—	—	—	—	02206 (4)	1	Adapter also available in pkg. of 1 00609 .

*Sealing Assembly replacement parts may be available. See section E for availability.

B-1