



S-20/20

Aluminum Body
 Titanium Buckets
 6-Place, 20 ml
 Tube Dimensions: 25 × 51 mm
 Swinging-Bucket
 Maximum Speed: 20,000 rpm
 K Factor: 369
 Rotor Mass: 6.7 kg

	RCF (g)	R (cm)
Max.	57,771	12.93
Min.	32,214	7.21

For use in RC-28S, RC-26 Plus, RC-24³, RC-5C Plus, RC-5C³, RC-5B Plus, RC-5B³, RC-5³, RC-2B⁴

Basic Catalog No. 11932	Complete Catalog No. 11947	Accessories
Includes: S-20/20 Rotor Vacuum Grease 65937 Lubricant 61556 Buckets, 20 ml, Set of 6 54293 O-ring, Viton®, 20 ml (6 supplied) 64667 Bucket Rack 52240 Rotor Extractor Tool 11946 Instruction Manual 11952	Includes Basic Catalog 11932 Plus: Polyallomer Tubes for 20 ml Bucket (8 pkgs, 25/pkg) 03288	Buckets ¹ , 17 ml, set of 6 52030 Buckets ¹ , 36 ml, set of 6 52035 Software Update Kit for RC-5C ² 12814

B-1

TUBES

< See page B-3 for Rotor to Centrifuge cross reference >

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	NOTES
NOMINAL	FILL					
20	20	Polyallomer Thin-Walled Tube	20,000	25 × 51	03288 (50)	Must be run full volume only.

⁽¹⁾Buckets are sold as a matched set and must always be used together. Interchanging of used buckets from one rotor body to another will invalidate warranty.

NOTE: Complete sets of buckets may be ordered separately for use on the S-20 rotor. Warranty becomes equal to that which remains on the rotor body.

⁽²⁾Contact your SORVALL® representative to determine your need for a Software Update Kit.

⁽³⁾Instrument no longer available.

All trademarks are the property of their respective owners.