

Thermo Scientific Dharmacon®

DharmaFECT® siRNA

Transfection Reagents

Product Description

- Panel of four transfection reagents optimized for efficient delivery of siRNA into mammalian cells.
- Provided as a liquid at a concentration of 1 mg/mL.
- Functionally verified to produce at least 75% transfection efficiency by assessing mRNA knockdown of positive control siRNAs.
- Confirmed to be free of microbial contaminants and ribonucleases

Shipping and Storage

- DharmaFECT siRNA Transfection Reagents are shipped on ice packs. If the ice melts during transit there is no risk to the product. Stability testing demonstrates no loss in product functionality after prolonged storage at warm temperatures (see the Stability Testing Product Bulletin for additional details).
- Upon receipts, store at 4°C. DO NOT FREEZE. DharmaFECT transfection reagents are labeled with a date of manufacture and are good for 18 months following that date when stored as recommended.

Related Products

- Dharmacon recommends including a positive and negative control, such as Dharmacon RNAi Control reagents, in every RNAi experiment. For more information, go to www.thermo.com/DharmaconRNAiControls

Accompanying Documents

- Basic siRNA Transfection Protocol

Publication Reference Guide

When referencing the use of DharmaFECT siRNA Transfection Reagents, please include the following information: DharmaFECT® (1, 2, 3, or 4) siRNA Transfection Reagent, Thermo Fisher Scientific, Lafayette, CO.

Contact Information

For technical questions regarding the use of siRNA reagents, please contact Dharmacon Products Technical Support at:

In America/Asia

800-235-9880

303-604-9499

dharmacon.lab@thermofisher.com

In Europe/Israel

00800 73724648

32-53-85-71-84

perbio.eurotech@thermofisher.com

In Other Countries

Please contact your appropriate distributor as listed on

www.thermo.com/dharmacondistributors

Product	Size	Catalog Number
DharmaFECT 1	0.2 mL	T-2001-01
	0.75 mL	T-2001-02
	1.5 mL	T-2001-03
	1.5 mL x 5 tubes	T-2001-04
DharmaFECT 2	0.2 mL	T-2002-01
	0.75 mL	T-2002-02
	1.5 mL	T-2002-03
	1.5 mL x 5 tubes	T-2002-04
DharmaFECT 3	0.2 mL	T-2003-01
	0.75 mL	T-2003-02
	1.5 mL	T-2003-03
	1.5 mL x 5 tubes	T-2003-04
DharmaFECT 4	0.2 mL	T-2004-01
	0.75 mL	T-2004-02
	1.5 mL	T-2004-03
	1.5 mL x 5 tubes	T-2004-04
DharmaFECT 1, 2, 3, and 4	0.2 mL each	T-2005-01
	0.75 mL each	T-2005-02
	1.5 mL each	T-2005-03

© 2008 Thermo Fisher Scientific Inc. All rights reserved. QuantiGene is a trademark of Bayer. alamarBlue is a trademark of Accumed International, Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Literature Code: 00089-06-F-02-U