

| CENTRIFUGES | | | | | | | | |
|-----------------------------|--------|-------------------------|---------------------|-------------------------------|-------------------------------------------------|--------------------------------------|-------------------------------------|----------------------------------------------------------------|
| MAXIMUM SPEED (rpm) | | | | | | | | |
| SUPERSPEED ROTORS | RC-28S | Evolution RC RC-26 Plus | RC-24 ³ | RC-5C Plus RC-5C ³ | RC-5B Plus RC-5 ³ RC-5B ³ | RC-2 ³ RC-2B ³ | SS-3 ³ SS-4 ³ | RC-3B, ³ RC-3B PLUS, RC-3C, ³ RC-3C PLUS |
| SE-12 | 21,000 | 26,000 | 24,000 | 21,000 | 21,000 | 21,000 | 16,000 | 6,000 ² |
| SA-300 | — | 25,000 ¹ | 24,000 | 20,500 ¹ | 20,500 | — | — | — |
| SS-34 | 20,000 | 20,500 | 20,500 | 20,000 | 20,000 | 20,000 | 17,000 | 6,000 ² |
| SM-24 | 20,000 | 20,500 | 20,500 | 20,000 | 20,000 | 20,000 | 16,000 | 6,000 ² |
| F-20/MICRO | — | 20,000 | 20,000 | 20,000 ¹ | 20,000 | 20,000 | — | — |
| SA-600 | 16,500 | 17,000 | 17,000 | 16,500 | 16,500 | 16,500 | 13,400 | 6,000 ² |
| SLA-600TC | — | 13,000 ¹ | 13,000 ¹ | 12,500 ¹ | 12,500 | 12,500 | — | — |
| SLA-1000 | — | 16,500 ¹ | 16,500 ¹ | 16,500 ¹ | 16,500 | 16,500 | — | — |
| SLA-1500 | — | 15,000 | 15,000 ¹ | 14,500 ¹ | 14,500 | 14,500 | — | — |
| SLA-3000 | — | 12,000 | 12,000 ¹ | 11,000 ¹ | 11,000 | 11,000 | — | — |
| GSA | 13,000 | 13,000 | 13,000 | 13,000 | 13,000 | 13,000 | 9,000 | — |
| GS-3 | 9,000 | 9,000 | 9,000 | 9,000 | 9,000 | 9,000 | — | — |
| SH-MT | 20,000 | 20,000 ⁶ | 20,000 | 20,000 ¹ | 20,000 | 20,000 | — | — |
| SH-80 | 20,000 | 20,000 | 20,000 | 20,000 ¹ | 20,000 | 20,000 | — | — |
| SH-3000/Buckets | — | 4,700 | 4,700 ¹ | 4,350 ¹ | 4,350 | — | — | — |
| SH-3000/Microplate Carriers | — | 4,700 | 4,700 ¹ | 4,350 ¹ | 4,350 | — | — | — |
| HB-6 | 13,000 | 13,000 | 13,000 | 13,000 ¹ | 13,000 | 13,000 | — | — |
| HB-4 ⁵ | 13,000 | 13,000 | 13,000 | 13,000 | 13,000 | 13,000 | 9,000 | 6,000 ² |
| HS-4 | 7,000 | 7,500 | 7,500 | 7,000 | 7,000 | 7,000 ⁴ | — | 6,000 ² |
| SV-288 | 20,000 | 20,000 | 20,000 | 20,000 | 20,000 | — | — | — |
| SV-80 | 19,000 | 19,000 | 19,000 | 19,000 | 19,000 | — | — | — |
| TZ-28 | 20,000 | 20,000 | 20,000 | 20,000 | 20,000 | — | — | — |
| TZ-28/GK | 19,000 | 19,000 | 19,000 | 19,000 | 19,000 | — | — | — |
| SS-34/KSB | 20,000 | 20,000 | 20,000 | 20,000 | 20,000 | 20,000 | — | — |

⁽¹⁾May require software update kit.

⁽²⁾Requires adapter **55222**.

⁽³⁾Instrument no longer available.

⁽⁴⁾For RC-2B only.

⁽⁵⁾Rotor no longer available. ⁽⁶⁾ For RC-26 plus only.

B-2