

**ROTORS**

**TUBES, BOTTLES**

(Number Represents Number of Tubes per Rotor)

NOMINAL CAPACITY (ml)	MATERIAL <sup>1</sup>	STYLE <sup>2</sup>	DIMENSIONS (O.D. x L mm)	CATALOG NUMBER <sup>3</sup>	SE-12	SA-300	SS-34	SM-24	F-20/MICRO	SA-600	SLA-1000	GSA, SLA-1500	GS-3, SLA-3000	SH-WT	SH-80	HB-6	HB-4 <sup>4</sup>	HS-4	SV-80	SV-288	SH-3000	SLA-600TC	
0.25	MICROTUBE	—	—	—					32					120									
0.4	MICROTUBE	—	—	—					32					120									
0.5	MICROTUBE	—	—	—					32					80									
1	CAB	T	7 x 50	03103		24	32			48						24	16	64					
1.5	MICROTUBE	—	—	—	12	12	16	24	32	24				60		12	8	32				132	
2	CAB	T	7 x 75	03235		24	32			48						24	16	64					
3	PY	T	10 x 75	03100		12	16	24		24	20	30	30			12	8	20				132	
4	PC	FL	11 x 75	03104		12	16			24	48	72	72			12	8	48				132	
4	PP	FL	11 x 75	03105		12	16			24	48	72	72			12	8	48				132	
4	—	SG	10 x 75	—							48	72	72					48				132	
5	PY	T	12 x 75	03102	12			24			48	72	72					96				132	
5	PA	TW	13 x 51	03127															16				
5	PA	UC	13 x 56	03898															16				
5	—	SG	12 x 75	—							48	72	72					96				132	
7	PP	FL	13 x 100	03121		6	8	24		12						6	4	16				88	
7	PC	FL	13 x 100	03120		6	8	24		12						6	4	16				88	
7	—	V	13 x 100	—																		88	
10	PA	TW	18 x 50	03968											8								
10	PA	TW,FL	13 x 100	03124		6	8	24		12						6	4	16				88	
10	PC	OR	16 x 83	03020		6	8			12						6	4	16				20	
10	PP	OR	16 x 83	03929		6	8			12						6	4	16				20	
10	PY	CN,G	18 x 100	—		6	8			12						6	4	16				20	
10	PY	T	13 x 100	—		6	8	24		12		30	30			6	4	48				88	
10	—	V	16 x 100	—																		72	
12	PC	FL	16 x 100	03115		6	8			12	28	42	42			6	4	28				72	
12	PP	FL	16 x 100	03116		6	8			12	28	42	42			6	4	28				72	
12	PY	T	18 x 102	03109		6	8			12						6	4	16				20	
14	PC	FL	18 x 75	03246	12						20											20	
14	PP	FL	18 x 75	03247	12						20											20	
14	PY	T	18 x 120	03111							20	30	30					20				20	
15	—	TC,CN	17 x 120	—							20	30	30					20				32	12
15	CO	T	18 x 102	03286		6	8	24		12						6	4	16				48	
15	—	UR	17 x 120	—																		32	
15	—	V	16 x 127	—																		72	
16	PC	FL	18 x 100	03243		6	8	24		12						6	4	16				48	
16	PP	FL	18 x 100	03244		6	8	24		12						6	4	16				48	
20	PA	TW	25 x 51	03288																	8		
30	CO	T	24 x 105	03163		6	8			12	12	18	18			6	4	16				12	
30	PC	OR	25 x 94	03719		6	8			12								16					
30	PP	OR	25 x 94	03718		6	8			12								16					
30	—	TC,FS	25 x 100	—																		20	

<sup>(1)</sup> MATERIAL CODES

CAB—Cellulose Acetate Butyrate  
 PY—Pyrex® Glass  
 PC—Polycarbonate  
 PP—Polypropylene  
 PS—Polysulfone  
 PA—Polyallomer  
 CO—Corex® Glass

<sup>(2)</sup> STYLE CODES

T—Thick-Walled  
 FL—Flanged  
 TW—Thin-Walled  
 UC—ULTRACRIMP®  
 DS—DRY-SPIN Bottle  
 V—Vacutainer  
 SG—Standard Glass  
 OR—Oak Ridge Bottle  
 CN—Conical  
 G—Graduated  
 TC—Tissue Culture  
 B—Bottle  
 UR—Urinalysis

<sup>(3)</sup>If Catalog Number is not supplied, tube/bottle is not available from SORVALL.

<sup>(4)</sup>Rotor no longer available.

All trademarks are the property of their respective owners.

**TO ORDER, CALL 1-800-522-SPIN (7746)**

**TUBES, BOTTLES** **ROTORS**  
(Number Represents Number of Tubes per Rotor)

NOMINAL CAPACITY (ml)	MATERIAL <sup>1</sup>	STYLE <sup>2</sup>	DIMENSIONS (O.D. × L mm)	CATALOG NUMBER <sup>3</sup>	SE-12	SA-300	SS-34	SM-24	F-20/MICRO	SA-600	SLA-1000	GSA/SLA-1500	GS-3/SLA-3000	SH-MT	SH-80	HB-6	HB-4 <sup>4</sup>	HS-4	SV-80	SV-288	SH-3000	SLA-600TC	
33	PA	UC	25 × 90	<b>03899</b>																8			
35	PA	CN,OR	29 × 105	—		6	8			12	4	6	6			6	4	16			20	12	
35	PC	CN,OR	29 × 105	—		6	8			12	4	6	6			6	4	16			20	12	
36	PA	TW	25 × 87	<b>03141</b>																	8		
50	—	TC,CN	30 × 115	—							4	6	6					12				20	12
50	—	TC,FS	30 × 115	—																			12
50	SS	T	28 × 101	<b>00579</b>		6	8			12						6	4	16				20	
50	SS	FL	28 × 101	<b>00517</b>		6	8			12								16				20	
50	PC	OR	28 × 103	<b>03480</b>		6	8			12						6	4	16				20	
50	PP	OR	28 × 103	<b>03530</b>		6	8			12						6	4	16				20	
50	PA	TW	28 × 104	<b>03139</b>		6	8			12						6	4	16				20	
50	PC	FL	29 × 102	<b>03146</b>		6	8			12						6	4	16				20	
50	PP	FL	29 × 102	<b>03147</b>		6	8			12						6	4	16				20	
50	PA	UC	29 × 108	<b>03528</b>		6	8			12													
100	PY	B	44 × 137	<b>03230</b>							4	6	6									4	
125	CO	B	53 × 118	<b>03231</b>							4	6	6					4				4	
150	CO	B	53 × 132	<b>03229</b>							4	6	6									4	
150	PC	T	45 × 132	<b>03143</b>							4	6	6									4	
200	SS	T	61 × 124	<b>00521</b>								6 <sup>5</sup>	6 <sup>6</sup>										
200	PP	CN	60 × 135	<b>03094</b>																		4	
250	PA	OR	60 × 120	<b>54477</b>						4	6	6						4				4	
250	PC	OR	61 × 124	<b>03939</b>						4	6	6						4				4	
250	PP	OR	61 × 124	<b>03937</b>						4	6	6						4				4	
250	PP	DS	61 × 146	<b>03069</b>						4	6	6										4	
250	PC	DS	61 × 146	<b>03098</b>						4	6	6										4	
250	—	TC, CN	60 × 177	—																		4	
285	SS	B	61 × 140	<b>00530</b>								6 <sup>5</sup>	6 <sup>6</sup>										
315	SS	B	61 × 153	<b>00522</b>								6 <sup>5</sup>	6 <sup>6</sup>										
500	PC	OR	70 × 160	<b>03944</b>									6										
500	PP	OR	70 × 160	<b>03943</b>									6										
500	PP	DS	70 × 190	<b>06829</b>									6										
500	PP	TC	96 × 146	—																		4	
750	—	B	98 × 135	—																		4	

**(1) MATERIAL CODES**

CAB—Cellulose Acetate Butyrate  
 PY—Pyrex® Glass  
 PC—Polycarbonate  
 SG—Standard Glass  
 PP—Polypropylene  
 PS—Polysulfone  
 PA—Polyallomer  
 CO—Corex® Glass

**(2) STYLE CODES**

T—Thick-Walled  
 FL—Flanged  
 TW—Thin-Walled  
 UC—ULTRACRIMP®  
 OR—Oak Ridge Bottle

CN—Conical  
 G—Graduated  
 TC—Tissue Culture  
 B—Bottle  
 DS—DRY-SPIN Bottle  
 FS—Free Standing Conical

<sup>(3)</sup>If Catalog Number is not supplied, tube/bottle is not available from SORVALL.

<sup>(4)</sup>Rotor no longer available.

<sup>(5)</sup>Bottles fit GSA rotor only.

<sup>(6)</sup>Bottles fit GS-3 rotor only.

All trademarks are the property of their respective owners.

**TO ORDER, CALL 1-800-522-SPIN (7746)**

**B-2**