



SA-800
 Aluminum Rotor
 8-Place, 100 ml
 Tube Dimensions: 38 × 105 mm
 Fixed Angle: 20°
 Maximum Speed: 20,500 rpm
 K Factor: 560
 Rotor Mass: 6.8 kg
 Autoclavable: 121°C

	RCF (g)	R (cm)
Max.	49,054	10.45
Min.	19,340	4.12

Catalog No. 76940

Includes:

- SA-800 Rotor
- Rotor Manual **76941**
- Rotor Extractor **76942**
- Cover O-ring (6) **20290483**
- Cover Locking Knob O-ring (6) **20290468**
- Lubricant **76003500**

TUBES, BOTTLES

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
94	85	Steel Tube	20,500	38 x 105	75002883 (1)	75001569 (1)	Sealing Cap	—		
94	70	Polyallomer, Thick Tube	20,500	38 x 105	76002872 (1)	75001568 (1)	Sealing Cap	—		
94	70	Polycarbonate Screw Cap Tube	20,500	38 x 105	75002829 (1)	75001540 (1)	—	—		
94	70	Polycarbonate Tube	20,500	38 x 105	76002810 (1)	75001568 (1)	Sealing Cap	—		
85	65	Polycarbonate Oak Ridge Tube	20,500	38 x 105	75002976 (1)	75002996 (1)	—	—		
50	28	Polypropylene Flanged Tube	20,500	29 x 102	03147 (25)	03268 (25)	Polypropylene, Snap-On	75003014	1	
50	28	Polycarbonate Flanged Tube	20,500	29 x 102	03146 (25)	03268 (25)	Polypropylene, Snap-On	75003014	1	
50	41	Polypropylene Oak Ridge Bottle	20,500	28 x 103	03530 (25)	03535 (25)	Polypropylene, Sealing	75003014	1	Nonsealing cap 03926 (25/pkg) available.
50	38	Polycarbonate Oak ridge Bottle	20,500	28 x 103	03480 (25)	03535 (25)	Polypropylene, Sealing	75003014	1	Nonsealing cap 03926 (25/pkg) available.
35	30	Polyallomer Nalgene Conical Oak Ridge Tube	5000	29 x 105	—	—	—	See Notes	1	Requires Nalgene adapter 3147-0050 .
35	30	Polycarbonate Nalgene Conical Oak Ridge Tube	5000	29 x 105	—	—	—	See Notes	1	Requires Nalgene adapter 3147-0050 .

(continued)

TUBES, BOTTLES

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
30	28	Polycarbonate Oak Ridge Bottle	20,500	25 x 94	03719 (25)	03614 (25)	Polypropylene, Sealing	00419 + 75003014 (1)	1	Nonsealing Cap 03925 (25/pkg) available.
30	28	Polypropylene Oak Ridge Bottle	20,500	25 x 94	03718 (25)	03614 (25)	Polypropylene, Sealing	00419 + 75003014 (1)	1	Nonsealing Cap 03925 (25/pkg) available.
16	12	Polyallomer Tube	20,500	18 x 100	75002951 (10)	75002957 (10)	—	00363 + 75003014	1	
16	12	Polycarbonate Flanged Tube	20,500	18 x 100	03243 (50)	03269 (50)	Polypropylene, Snap-on	00382 + 75003014 (1)	1	
15	10	Corex Tube	10,500	18 x 102	03286 (25)	—	—	00363 + 75003014 (1)	1	
12	6	Pyrex Tube	5000	18 x 102	03109 (50)	—	—	00363 + 75003014 (1)	1	
12	6	Polycarbonate Flanged Tube	20,500	16 x 100	03115 (50)	03266 (50)	Polypropylene, Snap-on	00402 + 75003014 (1)	1	
10		Pyrex Graduated Tube, Conical	3000	18 x 100	—	—	—	00367 (1)	1	
10	9	Polypropylene Oak Ridge Bottle	20,500	16 x 83	03929 (50)	03279 (25)	Polypropylene, Sealing	00425 + 75003014 (1)	1	Nonsealing Cap 03924 (25/pkg) available.
10	9	Polycarbonate Oak Ridge Bottle	20,500	16 x 83	03020 (25)	03279 (25)	Polypropylene, Sealing	00425 + 75003014 (1)	1	Nonsealing Cap 03924 (25/pkg) available.
10		Pyrex Tube	5000	13 x 100	—	—	—	00366 (1)	1	
7	4	Polycarbonate Flanged Tube	20,500	13 x 100	03120 (50)	03265 (50)	Polypropylene, Snap-on	00473 + 75003014 (1)	1	
4	2	Polycarbonate Flanged Tube	20,500	11 x 75	03104 (50)	03264 (50)	Polypropylene, Snap-on	00381 + 75003014 (1)	2	
3	1	Pyrex Tube	5000	10 x 75	03100 (50)	—	—	00364 + 75003014 (1)	2	
2	1	Cellulose Acetate Butyrate Tube	20,500	7 x 75	03235 (50)	—	—	00410 + 75003014 (1)	4	
1	1	Cellulose Acetate Butyrate Tube	20,500	7 x 50	03103 (50)	—	—	00408 + 75003014 (1)	4	