



**SA-512**  
 Aluminum Rotor  
 32-Place, 16 ml  
 Tube Dimensions: 18 × 100 mm  
 Fixed Angle: 23°  
 Maximum Speed: 19,500 rpm  
 K Factor: Outer 394, Inner 492  
 Rotor Mass: 7.3 kg  
 Autoclavable: 121°C

	RCF (g)	R (cm)
<b>Outer Row</b>		
Max.	48,930	11.52
Min.	27,056	6.37
<b>Inner Row</b>		
Max.	41,794	9.84
Min.	19,963	4.70

**Catalog No. 76950**

Includes:

- SA-512 Rotor
- Rotor Manual **76941**
- Rotor Extractor **76942**
- Cover O-ring (6) **20290484**
- Cover Locking Knob O-ring (6) **20290468**
- Lubricant **76003500**

**TUBES, BOTTLES**

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
16	12	Polyallomer Tube	19,500	18 x 100	<b>75002951</b> (10)	<b>75002957</b> (10)	Polypropylene Stopper	—	—	
16	12	Polycarbonate Flanged Tube	19,500	18 x 100	<b>03243</b> (50)	<b>03269</b> (50)	Polypropylene, Snap-on	—	—	
15	8	Corex Tube	10,500	18 x 102	<b>03286</b> (25)	—	—	<b>03667</b> (12)	1	
10		Pyrex Tube	—	13 x 100	—	—	—	<b>00373</b> (1)	1	
7	4	Polycarbonate Flanged Tube	19,500	13 x 100	<b>03120</b> (50)	<b>03265</b> (50)	Polypropylene, Snap-on	<b>00416</b> (1)	1	
5	3	Pyrex Tube	5000	12 x 75	<b>03102</b> (50)	—	—	<b>00376</b> (1)	1	
3	1.5	Pyrex Tube	5000	10 x 75	<b>03100</b> (50)	—	—	<b>00374</b> (1)	1	
1.5		Microtube	—	—	—	—	—	<b>00376</b> (1)	1	
1.5		Polypropylene Tube	19,500	11 x 38	<b>76001148</b> (500)	—	Incl. Screw Cap	<b>75002900</b>	1	