

SORVALL® SUPERSPEED CENTRIFUGES

SLC-6000 - SUPER-LITE® Composite Rotor



For use in: Evolution™ RC

SUPER-LITE® Composite Rotor
 6-Places, 1,000 ml
 Cavity Dimensions: 98 x 175
 Fixed Angle: 20°
 Maximum Speed: 8,500 rpm
 K Factor: 3,549
 Rotor Mass: 17 kg
 Autoclavable: Yes

	RCF (x g)	R (cm)
Max.	15,810	19.59
Min.	4,374	5.42

Catalog No. 76640	Catalog No. 76650	Catalog No. 74815
SLC-6000 Rotor Rotor Cover Assembly 78595 Cover O-ring 92939 Rotor Manual 76644 Rotor Extractor Tool 76645 Vacuum Grease 65937 Lubricant 61556 Polycarbonate Bottles (2 pkgs., 2/pkg.) Sealing Assembly (2 pkgs., 2/pkg.) 76960	SLC-6000 Rotor Rotor Cover Assembly 78595 Cover O-ring 92939 Rotor Manual 76644 Rotor Extractor Tool 76645 Vacuum Grease 65937 Lubricant 61556 Polypropylene Bottles (2 pkgs., 2/pkg.) 76968 Sealing Assemblies (2 pkgs., 2/pkg.) 76976	SLC-6000 Rotor Rotor Cover Assembly 78595 Cover O-ring 92939 Rotor Manual 76644 Rotor Extractor Tool 76645 Vacuum Grease 65937 Lubricant 61556 Wide-Mouth Polypropylene Bottle (2 pkgs., 2/pkg.) 74812 Sealing Assemblies (2 pkgs., 2/pkg.) 74813

Tubes and Adapters

CAPACITY (ML)		DESCRIPTION	MAX. SPEED (RPM)	TUBE DIMENSIONS (MM)	TUBE CAT. NO. (QTY/PKG)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT NO. (QTY/PKG)	DESCRIPTION	CAT. NO. (QTY/PKG)	NO. OF PLACES	
1,000	930	Polycarbonate Oak Ridge Bottles	8,500	98 x 175	76964 (2)	76960 (2)	Polyetherimide Sealing			
1,000	1000	Polypropylene Tubes	7,200	98 x 175	76968 (2)	76976 (2)	Polyetherimide Sealing			
1,000	900	Wide-Mouth Polypropylene Bottle	8,500	97 x 156	74812 (2)	74813 (2)	Polypropylene Sealing			
500	425	Polycarbonate Oak Ridge Bottles	8,500	70 x 160	03944 (6)	03280 (12)	Polypropylene Sealing	41003 (1)	1	Non-Sealing Cap available (12) (03928)
500	425	Polypropylene Oak Ridge Bottles	8,500	70 x 160	03943 (6)	03280 (12)	Polypropylene Sealing	41003 (1)	1	Non-Sealing Cap available (12) (03928)
500	480	Polypropylene DRY-SPIN™ Bottle	8,500	70 x 170	06829 (6)	Included with Bottles	Polypropylene Sealing	41003 (1)	1	
250	260	Polypropylene Oak Ridge Bottles	8,500	61 x 124	03937 (12)	03278 (12)	Polypropylene Sealing	41003 (1) & 41002 (1)	1	Non-Sealing Cap available (12) (03927)
250	240	Polycarbonate Oak Ridge Bottles	8,500	61 x 124	03939 (12)	03278 (12)	Polypropylene Sealing (12)	41003 (1) & 41002 (1)	1	Non-Sealing Cap available (12) (03927)

All trademarks are the property of their respective owners.

TO ORDER, CALL 1-800-522-7746 (SPIN)

SORVALL® SUPERSPEED CENTRIFUGES

Tubes and Adapters Continued for SLC-6000 Rotor

CAPACITY (ML)		DESCRIPTION	MAX. SPEED (RPM)	TUBE DIMENSIONS (MM)	TUBE CAT. NO. (QTY/PKG)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT NO. (QTY/PKG)	DESCRIPTION	CAT. NO. (QTY/PKG)	NO. OF PLACES	
250	260	Polycarbonate DRY-SPIN™ Bottle	8,500	61 x 146	03098 (12)	Included with Bottles	Polypropylene Sealing	41003 (1) & 41002 (1)	1	Replacement Caps (12) (06841)
250	250	Polypropylene DRY-SPIN™ Bottle	8,500	61 x 146	03069 (12)	Included with Bottles	Polypropylene Sealing	41003 (1) & 41002 (1)	1	Replacement O-ring sealing assembly (24) (03545)
50	41	Polypropylene Oak Ridge Bottles	8,500	28 x 103	03530 (25)	03535 (25)	Polypropylene Sealing	00894 (1)	5	Non-Sealing Cap available (12) (03926)
50	38	Polycarbonate Oak Ridge Bottles	8,500	28 x 103	03480 (25)	03535 (25)	Polypropylene Sealing	00894 (1)	5	Non-Sealing Cap available (12) (03926)

All trademarks are the property of their respective owners.

TO ORDER, CALL 1-800-522-7746 (SPIN)