

TUBES AND BOTTLES FOR SH-3000 (continued)

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
200	200	Nunc Tissue Culture Bottle, Conical	4,700	60 × 135	03094 (4)	Included with Bottle	Polypropylene Nonsealing	01092 (1) 03095 (1)	1 1	
150	144	Polycarbonate Thick-Walled Tube	4,700	45 × 132	03143 (25)	—	—	00458 (1) 01092 (1)	1 1	
150	120	Corex® Bottle	4,700	53 × 132	03229 (1)	Included with Bottle	Polypropylene Nonsealing	01092 (1) 00372 (1)	1 1	Bottle no longer available
125	100	Corex® Bottle	4,700	53 × 118	03231 (1)	Included with Bottle	Polypropylene Nonsealing	01092 (1) 00372 (1)	1 1	Bottle no longer available
100	100	Pyrex® Bottle	3,000	44 × 137	03230 (1)	Included with Bottle	Polypropylene Nonsealing	01092 (1) 00371 (1)	1 1	
50	—	Tissue Culture Tube, Conical	—	30 × 115	—	—	—	00436 (1) 00445 (1)	5 3	Use 00445 with sealing cover. Both adapters accommodate free standing and conical tissue culture tubes.
50	46	Polypropylene Flanged Tube	4,000	29 × 102	03147 (25)	03268 (25)	Polypropylene, Snap-on	00894 (1)	5	
50	48	Polycarbonate Flanged Tube	4,000	29 × 102	03146 (25)	03268 (25)	Polypropylene, Snap-on	00894 (1)	5	
50	56	Polyallomer Thin-Walled Tube	4,000	28 × 104	03139 (25)	—	—	00894 (1)	5	
50	41	Polypropylene Oak Ridge Bottle	4,000	28 × 103	03530 (25)	03535 (25)	Polypropylene, Sealing	00894 (1)	5	
50	38	Polycarbonate Oak Ridge Bottle	4,000	28 × 103	03480 (25)	03535 (25)	Polypropylene, Sealing	00894 (1)	5	
50	50	Stainless Steel Tube	4,000	28 × 101	00579 (1)	00518 (1)	Stainless Steel, Sealing	00894 (1)	5	Requires Wrench 01014 , for tube sealing.
50	50	Stainless Steel Flanged Tube	4,000	28 × 101	00517 (1)	00518 (1)	Stainless Steel, Sealing	00894 (1)	5	Requires Wrench 01014 , for tube sealing.
50	—	Flat-Bottomed Tube	—	29 × 100	—	—	—	00894 (1)	5	
35	—	Polyallomer Nalgene® Conical Oak Ridge Tube	—	29 × 105	—	—	—	00436 (1)	5	
35	—	Polycarbonate Nalgene® Conical Oak Ridge Tube	—	29 × 105	—	—	—	00436 (1)	5	
30	—	Flat-Bottomed Tube	—	25 × 100	—	—	—	00894 (1)	5	
30	31	Corex® Bottle	4,700	24 × 105	03163 (10)	—	—	00449 (1) 01092 (1)	3 1	Bottle no longer available
16	16	Polypropylene Flanged Tube	4,700	18 × 100	03244 (50)	03299 (50)	Polypropylene Nonsealing	00892 (1)	12	
16	17	Polycarbonate Flanged Tube	4,700	18 × 100	03243 (50)	03299 (50)	Polypropylene Nonsealing	00892 (1)	12	
15	—	Tissue Culture Tube, Conical	—	17 × 120	—	—	—	00447 (1)	8	
15	—	Urinalysis Tube, Conical	—	17 × 120	—	—	—	00447 (1)	8	
15	—	Vacutainer Tube	—	16 × 127	—	—	—	00385 (1)	18	
15	18	Corex® Tube	4,000	18 × 102	03286 (25)	—	—	00892 (1)	12	Rubber Pads 03667 (12/pkg), recommended.

*Sealing Assembly replacement parts may be available. See section E for availability.

(continued on next page)

All trademarks are the property of their respective owners.

TO ORDER, CALL 1-800-522-SPIN (7746)

B-2

TUBES AND BOTTLES FOR SH-3000 (continued)

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		ADAPTER REQ'D		NOTES
NOM-INAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty.Pkg)	NO. OF PLACES	
14	16	Pyrex® Tube	4,000	18 × 120	03111 (50)	—	—	00894 (1)	5	
14	12	Polypropylene Flanged Tube	4,000	18 × 75	03247 (50)	03269 (50)	Polypropylene, Snap-on	00894 (1)	5	5 plastic adapters 00382 required with 00894 .
14	12	Polycarbonate Flanged Tube	4,000	18 × 75	03246 (50)	03269 (50)	Polypropylene, Snap-on	00894 (1)	5	5 plastic adapters 00382 required with 00894 .
12	14	Pyrex® Tube	4,000	18 × 102	03109 (50)	—	—	00894 (1)	5	5 rubber adapters 00363 required with 00894 .
12	10	Polypropylene Flanged Tube	4,700	16 × 100	03116 (50)	03266 (50)	Polypropylene Nonsealing	00385 (1)	18	
12	11	Polycarbonate Flanged Tube	4,700	16 × 100	03115 (50)	03266 (50)	Polypropylene Nonsealing	00385 (1)	18	
10	—	Vacutainer Tube	—	16 × 100	—	—	—	00385 (1)	18	
10	—	Pyrex® Graduated Tube, Conical	—	18 × 100	—	—	—	00894 (1)	5	5 rubber adapters 00367 required with 00894 .
10	11	Polyallomer Thin-Walled Flanged Tube	4,700	13 × 100	03124 (50)	—	—	00384 (1)	22	
10	9.1	Polypropylene Oak Ridge Bottle	4,000	16 × 83	03929 (50)	03279 (25)	Polypropylene, Sealing	00894 (1) 00402 (5)	5 1	Nonsealing Cap 03924 (25/pkg) available.
10	9.7	Polycarbonate Oak Ridge Bottle	4,000	16 × 83	03020 (25)	03279 (25)	Polypropylene, Sealing	00894 (1) 00402 (5)	5 1	Nonsealing Cap 03924 (25/pkg) available.
10	—	Glass Tube	—	13 × 100	—	—	—	00384 (1)	22	
10	—	Pyrex® Tube	—	13 × 100	—	—	—	00384 (1)	22	
7	—	Vacutainer Tube	—	13 × 100	—	—	—	00384 (1)	22	
7	7	Polypropylene Flanged Tube	4,700	13 × 100	03121 (50)	03265 (50)	Polypropylene Nonsealing	00384 (1)	22	
7	7	Polycarbonate Flanged Tube	4,700	13 × 100	03120 (50)	03265 (50)	Polypropylene Nonsealing	00384 (1)	22	
5	—	Standard Glass Tube	—	12 × 75	—	—	—	00857 (1)	33	
5	—	Vacutainer® Tube	—	13 × 75	—	—	—	00857 (1)	33	Only 14 tubes per adapter due to size of plug.
5	—	Vacutainer® Tube	—	10.25 × 103	—	—	—	00857 (1)	33	
5	5	Pyrex® Tube	4,700	12 × 75	03102 (50)	—	—	00857 (1)	33	
4	—	Standard Glass Tube	—	10 × 75	—	—	—	00857 (1)	33	
4	—	Vacutainer® Tube	—	10.25 × 82	—	—	—	00857 (1)	33	
4	4	Polypropylene Flanged Tube	4,700	11 × 75	03105 (50)	03264 (50)	Polypropylene Nonsealing	00857 (1)	33	
4	4	Polycarbonate Flanged Tube	4,700	11 × 75	03104 (50)	03264 (50)	Polypropylene Nonsealing	00857 (1)	33	
3	3	Pyrex® Tube	4,700	10 × 75	03100 (50)	—	—	00857 (1)	33	
—	—	Microplate, 3 standard or 4 Deepwell	—	—	—	—	—	See Notes	3 1	Carrier 11779 is used in place of bucket.
1.5	—	Microtubes	—	—	—	—	—	00857 (1) 04871 (1)	33 33	Adapters 00857 and 04871 must be used together.

*Sealing Assembly replacement parts may be available. See section E for availability.

All trademarks are the property of their respective owners.

TO ORDER, CALL 1-800-522-SPIN (7746)

B-2