



SV-80

Aluminum Rotor
 16-Place, 5 ml
 Tube Dimensions: 13 x 51 mm
 Vertical: 0°
 Maximum Speed: 19,000 rpm
 K Factor: 98
 Rotor Mass: 6.9 kg

	RCF (g)	R (cm)
Max.	40,969	10.16
Min.	35,646	8.84

		RC-28S	RC-26 Plus	RC-24 ¹	RC-5C Plus RC-5C ¹	RC-5B Plus RC-5B ¹ RC-5 ¹	RC-2B ¹ RC-2 ¹	SS-3 ¹ SS-4 ¹	RC-3B ¹ RC-3B Plus RC-3C ¹ RC-3C Plus
SV-80	Max Speed (rpm)	19,000	19,000	19,000	19,000	19,000	n/a	n/a	n/a
	Max RCF (g)	40,969	40,969	40,969	40,969	40,969	n/a	n/a	n/a

B-2

Basic Catalog No. 09204	Complete Catalog No. 09200	Accessories
Includes: SV-80 Rotor Lubricant 61556 Rotor Extractor 01031 Instruction Manual 09107	Includes Basic Catalog 09204 Plus: ULTRACRIMP® Tubes, 5 ml, Polyallomer (4 pkgs, 50/pkg) 03898 ² Rotor Cap (16 pkgs, 1/pkg) 03092 Tube Racks (1 pkg, 2/pkg) 03768 ULTRACRIMP® Tube Removal Tool 12826	ULTRACRIMP® Sealing Tool and Crimp Gauge 03920 Crimp Gauge, Replacement 03919 ULTRACRIMP® Plugs and Caps, Extra (25 each/pkg) 03999 Tool Kit for ULTRACRIMP® 12770 . Includes: Rotor Holding Fixture 52540 , Torque Wrench 52580 , Socket 52882 Cap Wrench (for Polyallomer Thin-Walled Tubes) 52363 Tube Trimmer, optional 52387 Rotor Holding Fixture 52540 ³

TUBES

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	
5	5	Polyallomer ULTRACRIMP® Tube	19,000	13 x 56	03898 (50)	Included with Tubes	Plug and Aluminum Cap	Requires Rotor Cap 03092 , Tube Racks 03768 , and ULTRACRIMP® Sealing Tool 03920 , and Tool Kit 12770 .
5	4.4	Polyallomer Thin-Walled Tube	19,000	13 x 51	03127 (50)	09116 (1) 52522 (8)	Body Cap Plug	Requires Body Cap 09116 , Tube Plug 52522 , Rotor Holding Fixture 52540 , and Wrench 52363 .

⁽¹⁾Instrument no longer available.

⁽²⁾Tubes supplied with equal number of plugs and caps.

⁽³⁾Included in Tool Kit **12770** for ULTRACRIMP® tubes. Must be ordered separately for use with tubes **03127**.

*Sealing Assembly replacement parts may be available. See section E for availability.